

Expression and purification of CENP-T proteins

Expression

Day 1: Transform pET28a plasmid with CENP-T construct into ArcticExpress (DE3) *Escherichia coli* (Agilent Technologies, 230192). Grow culture in 150 mL LB/Kan media (LB + 50 µg/mL Kanamycin). Grow culture overnight at 37 °C / 240 rpm.

Day 2:

- 1) Split 150 mL of overnight culture into 2 x 1L of LB/Kan.
- 2) Incubate the cultures at 37 °C / 240 rpm.
- 3) Measure changes in OD600 over time.
- 4) Induction of protein growth:
 - a) When OD600 is about 0.3-0.4, transfer the culture to the cool shaker to 10 °C.
 - b) After 15-20 minutes, check OD600. When OD600 = 0.5-0.6 induce using 0.5 mM IPTG.
 - c) Incubate at 10 °C / 240 rpm for 22 hours.

Day 3:

- 1) Spin the cultures at 4,000 g / 4 °C / 30 minutes / (1L bottles, rotor Beckman JLA-9.1) and collect pellets.
- 2) Freeze the pellets in liquid nitrogen, and store them at -80 °C.

Purification

Protocol A. This version of the protocol was used to purify CENP-T^{6D}, CENP-T^{WT}, CENP-T^{2D}, CENP-T^{T11A}, CENP-TT^{T85A}, and CENP-T^{2D,short} proteins.

For purification of all CENP-T proteins, we use the fast protein liquid chromatography (FPLC) system AKTA Pure (GE Healthcare) at 4 °C.

1. Prepare buffers:

- 1) **Buffer A** (1L): 50 mM Tris-HCl pH 8.0, 300 mM NaCl, 5 mM β -mercaptoethanol. Use this buffer as a base for Lysis and Elution buffers.
- 2) **Lysis Buffer** (100 mL; prepare on the day of the experiment): Buffer A supplemented with 0.15 mg/mL lysozyme (Sigma-Aldrich, L6876), cOmplete EDTA-free protease inhibitors cocktail (Roche, 11873580001) and 1 mM phenylmethylsulfonyl fluoride (Sigma-Aldrich, P7626).
- 3) **Elution buffer** (250 mL): Buffer A supplemented with 500 mM imidazole; adjust pH after Imidazole dissolving.

2. Cell Lysis. All manipulations should be done on ice and/or in a cold room.

- 1) Thaw cell pellets on ice for 15 min (alternatively use fresh culture).
- 2) Resuspend cells in Lysis buffer (25 mL per 1 L of culture).
- 3) Sonicate on ice bath 30 sec ON, 30 sec OFF, 40% power (60 W), 2.5 min total time, $\frac{1}{2}$ in probe.
- 4) Spin the crude lysate by ultracentrifugation at 50,000 g for 30 min, 4 °C (Beckman 50.2Ti rotor). Collect supernatant.

3. Ni-NTA column. All manipulations should be done in a cold room.

- 1) Equilibrate 5 mL HisTrap column (Cytiva, 17524801) with 5 column volumes of Buffer A.
- 2) Transfer the cleared supernatant to the 100 mL super loop and attach it to AKTA.
- 3) Load the protein to the column at the flow rate of 2.5 mL/min.
- 4) Wash the column with ten volumes of Buffer A, five volumes of Buffer A supplemented with 25 mM imidazole, and five volumes of Buffer A with 50 mM imidazole at the flow rate of 2.5 mL/min.
- 5) Elute the protein with ten-column volumes of Buffer A containing the gradient of imidazole in the range 50-500 mM. The flow rate at this step is 2.5 mL/min. Collect 2 mL fractions.
- 6) Monitor the protein elution by absorbance at 280 nM.
- 7) Run the peak fractions on the SDS-PAGE to select those that contain clean CENP-T protein.

4. Size-exclusion chromatography. All manipulations should be done in a cold room.

- 1) Equilibrate to HiLoad Superdex 200 pg column (Cytiva, 28-9893-35) with 200 mL of Buffer A.

- 2) Combine the peak fractions containing CENP-T protein (~ 5 mL total), and concentrate them to 2-2.5 mL using 10 kDa Amicon Ultra Centrifugal Filter Units (EMD Millipore, UFC901008).
- 3) Centrifuged protein to remove aggregates at 30,000 g for 15 min, 4 °C.
- 4) Apply protein to the equilibrated HiLoad Superdex 200 pg column using a 5 mL loop and wash the loop with three loop volumes at the flow rate of 0.5 mL/min.
- 5) Elute the protein with Buffer A applied at a flow rate of 0.5 mL/min.
- 6) Monitor the protein elution by absorbance at 280 nM.
- 7) Run the peak fractions on the SDS-PAGE to select those that contain clean CENP-T protein. Based on the SDS-PAGE analysis, the fractions containing CENP-T proteins are in the 70-80 mL range.
- 8) Purified proteins were supplemented with 20% glycerol, aliquoted, flash-frozen in liquid nitrogen, and stored at -80 °C.

Protocol B. This version of the protocol was used to purify CENP-T^{site1}, CENP-T^{site2}, and CENP-T^{TR} proteins.

The protocol is similar to Protocol A with the following modifications:

1. Prepare 500 mL Buffer A with 0.5 mM EDTA.
2. Lysis Buffer had a different composition: Buffer A was supplemented with cOmplete protease inhibitors cocktail (Roche, 11697498001) and 0.5 mM EDTA.
3. Instead of step “**Ni-NTA column**” perform “**GST Sepharose purification**” described below.
4. All other steps of Protocol B are the same, as in Protocol A.

GST Sepharose purification. All manipulations should be done on ice and/or in a cold room.

- 1) Pack the gravity column with 3 mL of GST Sepharose 4B (Cytiva, 17-0756-01).
- 2) Equilibrate the column by washing with 20 mL of Buffer A with 0.5 mM EDTA.
- 3) Apply clear lysate using funnel extension to the funnel, let the column drop for 2 hours, and collect flowthrough.
- 4) After lysate has been applied to the column wash the column with 50 mL of wash buffer Buffer A with 0.5 mM EDTA.
- 5) Do an additional two washes by applying 50 mL of Buffer A with 0.5 mM EDTA each.
- 6) Close the outlet of the column with a cap.

- 7) Add to the column 2 mL of Buffer A with 0.5 mM EDTA supplemented with 200 μ L of 2 mg/mL of homemade TEV protease.
- 8) Close tightly the inlet of the column with a cap, additionally inlet and outlet of the column can be sealed with Parafilm.
- 9) Secure the column on the rotator and incubate overnight in a cold room (4 °C; 10-15 rpm).
- 10) Elute protein from the column by applying 2 mL of ice-cold Buffer A with 0.5 mM EDTA.
- 11) *If TEV cleavage reaction was not effective, close the cap of the column and add 2 more mL of Buffer A with 0.5 mM EDTA. Add 200 μ L of 2 mg/mL of homemade TEV protease. Incubate the column additional 1-1.5 hours at room temperature on the rotator (10-15 rpm). Elute protein by gravity followed by the addition of 2 mL of Buffer A with 0.5 mM EDTA.*