

**University of Pennsylvania
Perelman School of Medicine
High-Throughput Screening Core**

Sara Cherry, Ph. D.

David C. Schultz, Ph. D.

dschultz@mail.med.upenn.edu

(215) 573-9641

67 John Morgan Building



Mission

- **Provide the PSOM community with HTS resources.**
 - To educate and assist with HTS assay development, optimization, miniaturization, and validation
 - To provide laboratory robotics infrastructure and technically trained staff for HTS
 - To provide libraries of small molecule and genetic tools for HTS
 - To facilitate small-scale screens from user-defined gene-sets
- **Develop novel technology**
- **Seed collaborative research programs and grants.**
- **Educate the SOM on utility and uses of HTS**

SOM Screening Core

- › **Libraries**
 - › Maintenance
 - › Distribution
- › **Liquid handling**
 - › Janus MDT/Verispan 8-tip
 - › Bulk reagent dispensors
 - › ELx405 microplate washer
- › **Assay Detection**
 - › EnVision multi-mode microplate reader
 - › ImageXpress Micro
 - › FLIPR screening system
- › **BSL2 Tissue Culture capabilities**
 - › Hood
 - › Incubators
- › **Informatics**
 - › Automated image analysis
 - › Statistical analysis
 - › Validation analysis

Automation: Screening



Automated pippeting station




384w plates spotted with cDNAs or siRNAs





Dispenser



HT Multimode Reader




"Hit"



Automated Microscopy

Slide courtesy of J.Hogenesch

Cell-based functional screening approaches



SOM Screening Core Chemical Libraries

Bioactives, FDA approved, and FDA-like compounds

- SelleckChem Bioactives (~2100)
 - Kinase Inhibitors, Epigenetic Inhibitors, Cancer compounds, GPCR/Ion Channel, Metabolism, Microbiology, FDA approved/FDA-like
- LOPAC (1280): Library of Pharmacologically Active Compounds

Natural Products

- Microsource Purified Natural Products (800)

Diversity set

Chembridge (32,000 total)

20,000 from Core set and 12,000 from Express Pick set


ChemDiv (12,000 total)

SMART library

Chemical Libraries:

SelleckChem Bioactive Library (2083 cmpds)


•Kinase Inhibitors, Epigenetic Inhibitors, Cancer compounds, GPCR/Ion Channel, Metabolism, Microbiology, FDA approved/FDA-like



<http://www.selleckchem.com/screening/chemical-library.html>

Chemical Libraries: LOPAC¹²⁸⁰: Library of Pharmacologically Active Compounds

1,280 FDA approved, marketed drugs with annotated biological activities, predictable activities and proven scaffolds directed against a wide range of drug targets







<http://www.sigmaaldrich.com/catalog/product/sigma/lo1280?lang=en®ion=US>





Diversity Library: cherry picked 44K from...

- Started with ~200,000 compounds
- Used modified Lipinski parameters to filter the set.
- Removed Reactive groups (e.g. Michael acceptors) and compounds with undesirable functionalities.
- Duplicates were removed.
- Lead Finder Clustering and MACCS fingerprints were generated.
- PAINS1, PAINS2 and PAINS3 substructure filters were applied.
- REOS filters were applied.
- Property Based Selection of 50K cpds was performed
- Vendor re-supply validated

44K Diversity Set Characteristics



44K Diversity Set Characteristics-II



SOM Screening Core Genetic Libraries

siRNA

- human genome-wide, human drugable genome, human GO categories
- user-defined human and mouse

Non-coding RNAs

- lncRNAs (human)
- miRNA mimics/antagonists (human V20)

Human TRC 2.0 and Mouse TRC1.0 Lentivirus shRNA library

- Screening pools: GO categories; user-defined sets
- Order groups/individuals

MGC cDNA collection (CMV-driven)

- 18,000 full length, sequenced, mouse and human (arrayed);
- user-defined sets
- Order groups/individuals

Genetic Libraries: Ambion Silencer Select siRNA & miRNA mimics/inh.

A. Human genome Silencer Select siRNA

	Genes	siRNAs	siRNAs/GT	format	# plates*
Druggable genome [^]	9032	27093	3	pooled	26
Druggable genome Ext.	1383	4149	3	pooled	4
Rest of Genome	11170	33510	3	pooled	34

[^] Further organized into GO categories (e.g. kinases, NHRs, GPCRs)

B. Human non-coding RNA siRNA (lnc) & miRNA mimics/inhibitors

	Targets	siRNAs	siRNAs/Tar	format	# plates*
long non-coding RNAs (lnc)	2220	6660	3	pooled	8
miRNA mimics (human V20)	2555	2555	1	individual	8
miRNA inhibitors (human V20)	2555	2555	1	individual	8

* All assay plates are pre-spotted with siRNAs/miRNAs, including controls

Control siRNAs

Control	siRNA_ID	Type	Library
Neg1	S813	Negative	Genome, Inc
Neg2	S814	Negative	Genome, Inc
GAPDH	S815	Positive mammalian	Genome, Inc
GFP	s229097	Positive Non-mammalian	genome
Luciferase	s229095	Positive Non-mammalian	genome
Kif11	S7903	Positive Mammalian (death)	Genome, Inc, mir
MALAT1	s239370	Positive Mammalian	Lnc

* Controls are pre-spotted in all assay plates

Custom siRNA Libraries

- Custom siRNA libraries are provided by Ambion/Life technologies
- Minimum order of 20 siRNAs **REQUIRED**
- 0.1 nmol of Silencer Select siRNA (human)
- 1 nmol of Silencer siRNAs (human and mouse)
- siRNAs are received in microtiter plates (96 or 384)
- **All orders are submitted through the HTSC**

Custom Library Ordering

- **Complete the order form**
 - Customer information (26 digit Penn fund no) (**REQUIRED**)
 - NCBI Gene ID, Gene Symbol, # of siRNAs/ gene target, species (**REQUIRED**)
 - Email to dschultz@upenn.edu

Product Line (select from dropdown list)	siRNA or miRNA ID	# for same target*	Catalog Number	Entrez Gene ID**	Gene Symbol
Silencer® Select siRNA		3	4427030	23468	CBX5
Silencer® Select siRNA		3	4427030	1660	DHX9
Silencer® Select siRNA		3	4427030	3069	HDLBP
Silencer® Select siRNA		3	4427030	27316	RBMX, RBMXL1
Silencer® Select siRNA		3	4427030	23435	TARDBP
Silencer® Select siRNA		3	4427030	5725	PTBP1

Custom Library Ordering (cont)

➤ Costs

- **HTS/Omics validation set**

100 gene targets x 3 siRNAs is ~\$4,500

- **Custom library (Epigenetics/metabolism)**



500 gene targets x 3 siRNAs is ~\$12,000

➤ Turnaround

Order size	Order to deliver
20-250 siRNAs	5-15 business days
251-500 siRNAs	16-25 business days
501-1000 siRNAs	20-30 business days
1001-2500 siRNAs	25-30 business days

Genetic libraries: TRC human/mouse LV-shRNA library

		Gene Targets	clones	Validated	% val
Human	complete	20,018	129,695	43,470	34
Mouse	V1.0	15,960	77,819	20,564	26
Totals		35,978	207,514	64,034	



<http://www.sigmaaldrich.com/life-science/functional-genomics-and-rnai/shrna/library-information.html>

shRNA library services

➤ **Distribution of Individual clones**

- \$100.00 per Gene Target (5+ clones)
- Investigator should sequence confirm shRNA sequence
- Investigator responsible for validating knock-down

➤ **Preparation of Custom Gene Sets**


- Cherry-pick and array clones of interest into new plates
 - GO categories (e.g. Epigenetic targets)
 - User-defined sets (e.g. 'omics validation)
- Glycerol stock plate + pooled plasmid DNA
- 150 gene Targets/625 clones: \$2000
- 600 gene Targets/3000 clones: \$5000

Genetic libraries: MGC cDNA collection

- Obtained by J.Hogenesch
- 18,000 full length, sequenced, mouse and human cDNAs pre-cloned into pCMV-SPORT6
 - Insert is fully sequenced and guaranteed to match corresponding BC Accession Number
 - Expression-ready vectors eliminate additional cloning steps
 - Robust CMV promoters drive cDNA expression
 - Gateway sites flanking coding sequence allow for additional flexibility
 - Stored as bacterial glycerol stocks, arrayed in 96 well microtiter plates

Screen:

- Complete library
- user-defined sets (e.g. Interferon Stimulated Genes-ISGs)
- Order individual clones or groups
 - \$50.00/clone



MCS1: KpnI, AgeI, EcoRI, SmaI, EcoRV, and SalI

MCS2: NotI, XbaI, XhoI, ApaI, HindIII, and MluI

Automation: Screening



Automated pippeting station




384w plates spotted with cDNAs or siRNAs




Dispenser





HT Multimode Reader



*
* * *
* * *



"Hit"



Automated Microscopy

Slide courtesy of J.Hogenesch

Envision Xcite: multi-mode microplate reader

A. Measurement Technologies

- absorbance,
- fluorescence intensity (FI),
- fluorescence polarization (FP),
- time-resolved fluorescence (HTRF)
- ultra-luminescence
- AlphaScreen (Amplified Luminescence Proximity Homogeneous Assay)

B. Assays

Anything that requires a plate reader!

- Luciferase Reporter Gene Assays (RGAs)
- Viability assays (e.g. Cell-titer Glo)
- Enzymatic assays (FI)
- Protein: nucleic acid interaction (FP, HTRF, alphascreen)
- Protein: protein interactions (FP, HTRF, alphascreen)
- AlphaLISA




ImageXpress Micro XLS

A. Enables automated acquisition of multi-channel fluorescence images in microtiter plates

- 5 colors

B. Phenotypic assays upto 40X

- Physical properties of cells (size, shape, etc.)
- Puncta formation (autophagy, lipid droplets, etc)
- Activity (migration, invasion, neurite outgrowth)
- Expression and localization of native proteins
- Measurement of fluorescent labels (antibody stains, EdU/BrdU, Calcein, PI uptake, mitochondria, etc.)
- FRET




C. Image analysis: MetaXpress


- Easy analysis pipelines in-house
- Computational analysis of images enables measurement of multiple cellular properties of interest at the object/cell level and/or population level.

High-content Screening Assays


Neurite Outgrowth




Migration/scratch assay



3D tumor spheroids (Invasion/viability)



Lipid droplets/puncta



FLIPR Tetra: Fluorescence Imaging Plate Reader

A. Measurement Technology:

- Integrates liquid handling with rapid whole plate imaging of fluorescence and luminescence
- Can read at sub-second intervals, which enables the kinetics of the response to be captured
- Integrated pipettor enables successive liquid additions, providing an opportunity to detect agonists, antagonists, and allosteric modulators all in one assay


B. Assays:

1. Intracellular Ca²⁺ flux: cell permeable Ca²⁺ sensitive fluorescent dye (e.g. cardiotoxicity);
2. Fatty acid or neurotransmitter uptake of fluorescence dyes.
3. Membrane potential: lipophilic, anionic, fluorescent dye that partitions across the cytoplasmic membrane of live cells based membrane potential.




FLIPR Assays

Monitoring Intracellular Ca²⁺ Mobilization



The Assay Development Process



Assay Design and Development

- Is the nature of the response clearly defined?
- Is the response dependent only on the activity of the compound being tested or is it conditioned by another stimulus?
- What is the duration of the response?
- What sort of 'secondary assays' exist to confirm activity and determine selectivity of newly identified probes?
- What structural classes of known actives exist and what are the limitations to their use?

Goal: To design an assay that is robust and sensitive

What biology and assays are readily screened?

- Enzymatic
- Protein: protein interactions
- Protein: nucleic acid interactions
- Reporter Gene Assay
- Sensitivity to drugs
- Simple design and robust signals...



Biochemical Assays

- **Enzymatic assays**
 - Luminescence (e.g. ADP-Glo)
 - Fluorescence (e.g. QFRET, AmplexRed)
 - Transcreener
- **Protein-protein interactions**
 - ALPHAscreen
 - HTRF
 - ELISA (e.g. luminescence)
- **Protein-nucleic acid interactions**
 - ALPHAscreen
 - HTRF
- **Ligand binding**
 - Differential Scanning Fluorimetry/Thermal Stability
 - SPR

Profiling and Sensitivities

- Phenotypic profiling of cell lines
 - FDA and FDA-likes
 - Annotated gene family (e.g. kinome)
 - Synthetic lethality screens
 - Synergy studies (combinations gene-gene; gene-drug; drug-drug)
- Differentiated iPSC from individuals
- Across tumors (e.g. melanomas)
- Of a particular patient tumor line (to define responsiveness)

Pharmacological Response of Genetically Defined Tumor Cells



Cell-based Assays

- **Anything you can read in a plate reader**
 - Reporter Gene Assays (e.g. Luciferase)
 - Survival (e.g. Cell-titer Glo, Alomar Blue, MTS)
 - Signaling (e.g. alphascreen, alphaLISA)
- **Microscopy assays**
 - EdU/BrDU incorporation
 - Autophagy (e.g. LC3-GFP/RFP)
 - Migration assays (e.g. scratch assay)
 - Neurite Outgrowth
 - Infection
 - Nuclear/cytoplasmic shuttling
 - Antigen localization/staining
- **Signaling** (e.g. Ca^{2+} , membrane potential)

Cell types

- **Chemical screens- all cells amenable**
- **Genetic screens**
 - **Transformed cells**
 - **Many primary cells**
 - **Fibroblasts, endothelial, epithelial, macrophages, etc**
 - **Lymphocytes challenging**
 - **Stem cells**
 - **Mixed cultures**
 - **Perturb one cell type and assay another**
 - **e.g. siRNA in fibroblast but read tumor cell biology; siRNA in macrophage but read T cell biology**

What scale is useful to you?

- Do you have a defined set of genes or pathways you are interested in?
 - Functional study of ‘OMICS data’
 - RNAseq
 - CHIPseq
 - GWAS
 - Exome
 - Interested in Kinases or another ‘category’

Services

- **Consultations** (per hour)
 - Assay development
 - Assay optimization
 - Assay validation
 - Grant submissions
- **Equipment usage** (per hour)
 - With help
 - Without help
- **Small scale screens**
 - User-defined (siRNAs, shRNAs, cDNAs, chemicals)
 - Library plates (e.g. kinome)
- **Large scale screens**
 - Library plates
- **Data Analysis**
 - Normalization, annotation
 - HCA analysis sequence dev.
 - Screen reports
- **Reagents**
 - Transfection
 - Plastics
 - Tips
 - Plates
- **siRNAs, shRNAs, cDNAs**
 - User defined sets
 - Individual clones

Funding Opportunities

- NIH

- PAR-13-364 Development of Assays for High-Throughput screening for use in Probe and Pre-therapeutic Discovery (R01)
- PAR-14-283/PAR-14-284, High-Throughput Screening (HTS) to Discover Chemical Probes (R21/R01)
- PAR-14-279, Discovery of *in vivo* Chemical Probes (R01)
- PAR-13-049/PAR-13-048, Drug Discovery for Nervous System Disorders (R01/R21)
- PAR-13-007, Early-stage Pharmacological Validation of novel Targets and Accompanying Pre-therapeutic Leads for Diseases of Interest to the NIDDK (R01)
- PAR-14-006, Seeding Collaborations for Translational Research to Discover and Develop New Therapies for Diseases and Conditions within NIDDK's Mission (R01)
- PAR-15-056, Building on High Impact Neurobiology Through Assay Development: Advancing Tools for Therapeutic Development (R01)
- PAR-15-070, Innovation Grants to Nurture Initial Translational Efforts (IGNITE): Assay Development and Therapeutic Agent Identification and Characterization to Support Therapeutic Discovery (R21/R33)
- PAR-15-071, Innovation Grants to Nurture Initial Translational Efforts (IGNITE): Pharmacodynamics and In vivo Efficacy Studies for Small Molecules and Biologics/Biotechnology Products (R21/R33)

Funding Opportunities (part II)

- **NIH** (continued)
 - PAR-13-267, Novel NeuroAIDS Therapeutics: Integrated Preclinical/Clinical Program (P01)
 - PAR-15-041, Targeting Persistent HIV Reservoirs (TaPHIR) (R21/R33)
- **NCAT/TRND opportunities**
- **Foundations** (e.g. Welcome Trust, Melanoma Research Foundation, Leukemia/Lymphoma Society, Gates, Cystic Fibrosis)
- **Commercial** (e.g. Bayer Grants4targets, Astrazeneca Openinnovation)
- **Institute/Center/Program Pilot project funds**
 - Institute for Immunology
 - Center for Orphan Disease Research
 - Institute for Regenerative Medicine
 - Epigenetics Program

Input from the Community

We are expanding based on SOM needs....

- FLIPR from Physiology
- cDNA library from J. Hogenesch (Pharm)
- Mouse shRNA library
 - Contributions: CHOP, Cancer Biology and IRM
- Small molecule libraries
 - Diversity library contributions: CDB, Biochemistry and Microbiology
 - LOPAC from J. Hogenesch (Pharm)

What libraries or functionality would be useful to you?

HOW TO GET STARTED?

- Contact David Schultz at dschultz@mail.med.upenn.edu or 215-573-9641 for an initial consultation
 - Define the project
 - Determine if the facility has relevant expertise/technology to pursue the project
 - Develop a management plan
 - Set expectations
 - Get started!