

Digitonin treatment

- to disrupt the plasma membrane and wash out cytoplasmic, "unbound" proteins

Reference EMBO J., 20, 272-284 (2001)
"To remove cytosolic staining, cells were incubated in digitonin buffer (25 mM HEPES-KOH, pH 7, 125 mM KOAc, 2.5 mM Mg(OAc)₂, 0.004% digitonin) at room temperature for 2 minutes prior to fixation." 0.004%=40 µg/mL

Comment According to the result of HeLa cell expressing C-terminal fragment of mtGolgin-1, which localizes to TGN, I think 20 µg/mL digitonin was enough to see localization by immunofluorescence microscopy. Therefore you may reduce the concentration of digitonin to carry out this experiment.

Reagents

Digitonin stock solution

20 mg/mL digitonin in DMSO.
Stock at -20 C.

1 M HEPES-KOH, pH 7.0 (HEPES : FW 238.3)

Dissolve HEPES in 80% volume of water, adjust pH to 7.0 with KOH, and then volume up to 100% volume.

2 M KOAc (FW 98.14)

1 M Mg(OAc)₂ (FW 214.46)

KHM buffer

1 M HEPES-KOH, pH 7.0	x40 dilution
2 M KOAc	x16 dilution
1 M Mg(OAc) ₂	x400 dilution

Digitonin in KHM buffer

Dilute digitonin stock solution to a final concentration of 20~40 µg/mL with KHM buffer.

Protocol

- 0) Prepare cells on coverslip.
- 1) Wash cells with KHM buffer.
- 2) Add 20~40 µg/mL digitonin in KHM buffer and incubate at room temperature for 2 minutes. Remove the solution.
- 3) Fix the cells : add 2% formaldehyde and leave at room temperature for 30 minutes. Wash the cells twice with PBS for 10 minutes.