

Metabolic Labeling (a la Mickey)

Special Stuff you need to have long before you start:

Dialyzed FBS: Serum dialyzed exhaustively against hepes-buffered-saline (10mM HEPES pH 7.4/ 150mM NaCl) or PBS, and then sterile filtered, aliquoted, and stored at -20°C

1500X methionine/cysteine: 0.3M methionine and 0.3M cysteine in ddH₂O, sterile filtered and stored at -20°C

Met-/Cys- medium: DMEM without methionine or cysteine (from Gibco/BRL or other commercial source)

Tran35S-Label or other source of metabolic label

Reagents to prepare with the above shit and other stuff:

Met-/Cys- Complete Medium

DMEM without methionine or cysteine

3-5% dialyzed FBS

glutamine

pen/strep

10-25mM HEPES pH 7.4

Chase medium, warm and ice cold

Normal, complete medium with serum

15X methionine/cysteine (add 1/100th volume of 1500X met/cys)

10mM HEPES

Ice cold PBS

Hood space and a 37°C water bath

Procedure:

1. If using adherent cells, trypsinize cells from dish or flask in the usual manner;
 - resuspend in ~10ml of complete medium (with FBS), transfer to a centrifuge tube, and spin out cells at 1000-1500 RPM in the table-top centrifuge for 5-7'
 - if using non-adherent cells, simply spin them out: usually $5-10 \times 10^6$

2. Remove all traces of medium
 - resuspend cells in complete met-/cys- medium to $\sim 5 \times 10^6$ cells/ml (this means about 1 10cm dish in ~ 2 ml total).
 - incubate at 37°C for 15-30' for transient HeLa or COS cells (high expressors of whatever you are looking for); 30-60' for stables or other cell lines (low expressors of whatever you are looking for).

3. In the meantime, prepare the labeling medium:
 - for pulse / chases, I prefer to use $\sim 2\text{mCi/ml}$ of $\text{Tran}^{35\text{S}}$ -Label or other $^{35\text{S}}$ -methionine source; for long-term labeling (2 hrs or more), use 0.25 0.5mCi/ml .For a 10cm dish-ful of cells, 0.5-1ml is fine. (For these, you could actually just add the $^{35\text{S}}$ -met straight to the pre-incubating cells). For non-adherent cells, shoot for $5-10 \times 10^6$ per ml.
 - keep medium at 37°C until ready to add

4. At the end of the pre-incubation, spin out cells.
 - remove medium
 - resuspend in labeling medium, and incubate at 37°C for however long you want (a typical pulse is 30'); shake tube every now and then during the labeling
 - at this point, it may be convenient to aliquot cells into pulse and chase samples

5. At the end of the labeling:
 - cool cells on ice
 - spin cells out at 4°C
 - resuspend chase samples in warm chase medium (1ml or more is convenient, usually); incubate at 37°C for however long you want.
 - it may be convenient to separate multiple chase samples either at the beginning, i.e. here, or later on, depending on your preference and/or space in the water bath
 - resuspend pulse sample in ice cold chase medium as a 1st wash

6. At the end of the chase (and the pulse sample);
 - spin out cells at 4°C
 - remove medium; save if necessary
 - resuspend cells in a reasonable volume of ice cold PBS or Tris buffered saline, and spin out again
 - discard PBS wash and freeze cell pellets on dry ice for storage, or go ahead and lyse cells directly as desired for other procedures