

Penn Center for Musculoskeletal Disorders (P30-AR069619)

PCMD MicroCT Imaging Core Learning Lunch Series

Introducing our NEW μ CT45 scanner featuring a 20-sample holder carousel

X. Sherry Liu, Ph.D
Wei-Ju Tseng, M.S
Yilu Zhou, Ph.D
Rebecca Chung, Ph.D

Dec 4th, 2019

McKay Orthopaedic Research Laboratory
Department of Orthopaedic Surgery
Perelman School of Medicine
University of Pennsylvania
Philadelphia, PA



Penn


McKay Orthopaedic Research Laboratory



Personnel



X. Sherry Liu, Ph.D
Director



Yilu Zhou, Ph.D
Wei-Ju Tseng, MS
Manager



Rebecca Chung, Ph.D
Coordinator

Our contact email:
pcmd.microct@gmail.com

Available Scanners at PCMD μ CT Imaging Core


- Specimen μ CT

- μ CT 35
- μ CT 50
- μ CT 45 *new!*




- *In Vivo* μ CT

- vivaCT 40
- vivaCT 75
- vivaCT 80 *new!*



- Clinical μ CT

- XtremeCT II



Which Scanner Should I Use for My Study?

- μ CT 35 and μ CT 45 (McKay Lab)
 - Small specimen scans (under diameter of 35-45 μ m)
 - Mouse bone microstructure phenotyping **must** use μ CT 35/45
- vivaCT 40 (McKay Lab)
 - Rodent study (rats and mice) requiring in vivo scans (IACUC approval required)
 - When μ CT 35/45 is fully occupied
- vivaCT 80 (McKay Lab)
 - Rodent study (rats and mice) requiring in vivo scans (IACUC approval required)
 - In vivo study of rat vertebrae and skull must use vivaCT 80 instead of vivaCT 40
 - Large specimen scans (above diameter of 50 μ m)
- vivaCT 75 and μ CT 50 (VA Hospital)
 - Only accessible to investigators with VA affiliations
 - μ CT 50: studies requiring high resolution characterization
 - vivaCT 75: large specimen scans (above diameter of 50 μ m)
- Extreme CT II (CHOP Nutrition and Growth Lab)
 - Clinical studies (IRB approval required)
 - Large specimen scans (above diameter of 80 μ m)
- **For consistent results, please use the same model of scanner for all samples/animals from the same study**



PCMD *ex vivo* μ CT Scanners

	μ CT 35	μ CT 45	μ CT 50
Use	Specimen	Specimen	Specimen
X-Ray Source	30 - 70 kVp	45, 55, 70, 90 kVp	30 - 90 kVp
Max Scan Size	37.9 x 120 mm (\varnothing x L)	50 x 120 mm (\varnothing x L)	50 x 120 mm (\varnothing x L)
Max Specimen Size	75.8 x 140 mm (\varnothing x L)	90 x 120 mm (\varnothing x L)	100 x 160 mm (\varnothing x L)
Best image voxel size	3.5 μ m (\varnothing : 7 mm)	3.0 μ m (\varnothing : 9 mm)	1.5 μ m (\varnothing : 3 mm)
Location	McKay Lab	McKay Lab	VA Hospital

μ CT 35


Time cost per stack: ~23 min

Tube Diameter	Max Resolution	Slice# per Stack	Max length per stack
7 mm	3.5 μ m	235	0.81 mm
11.5 mm	6 μ m	233	1.40 mm
20 mm	10 μ m	232	2.32 mm
30 mm	15 μ m	231	3.47 mm
37 mm	18.5 μ m	231	4.27 mm


μ CT 45

Time cost per stack: ~48 min

Tube Diameter	Max Resolution	Slice# per Stack	Max length per stack
9 mm	3 μ m	1032	3.10 mm
14 mm	4.5 μ m	1032	4.64 mm
24 mm	7.4 μ m	1032	7.64 mm
34 mm	10.4 μ m	1032	10.73 mm
48 mm	14.6 μ m	1032	15.07 mm



Adjust stacks



μCT 35

Time cost per stack: ~23 min

Tube Diameter	Max Resolution	Slice# per Stack	Max length per stack
7 mm	3.5 μm	235	0.81 mm
11.5 mm	6 μm	233	1.40 mm
20 mm	10 μm	232	2.32 mm
30 mm	15 μm	231	3.47 mm
37 mm	18.5 μm	231	4.27 mm

μCT 45

Time cost per stack: ~48 min

Tube Diameter	Max Resolution	Slice# per Stack	Max length per stack
9 mm	3 μm	1032	3.10 mm
14 mm	4.5 μm	1032	4.64 mm
24 mm	7.4 μm	1032	7.64 mm
34 mm	10.4 μm	1032	10.73 mm
48 mm	14.6 μm	1032	15.07 mm

μCT 35: 3.5 μm for 1032 slices: 115 min

μCT 45: 3 μm for 1032 slices: 48 min

μCT 35: 10 μm for 1032 slices: 115 min

μCT 45: 10.4 μm for 1032 slices: 48 min

Resolution	You can use these tubes
3 μm	9mm
4.5 μm	9mm, 14mm
7.4 μm	9mm, 14mm, 24mm
10.4 μm	9mm, 14mm, 24mm, 34mm
14.6 μm	9mm, 14mm, 24mm, 34mm, 48mm


μ CT 35

6.0 μ m



μ CT 45


7.4 μ m



Penn


μ CT 35

6.0 μ m




μ CT 45

7.4 μ m





Phenotyping

- Mouse trabecular bone (long bone)
 - Distal femur, proximal tibia
 - μ CT35:
 - 6 μ m, 1-2 stacks, ~210-420 slices (1.4-2.8mm); Scan time: 24-48 mins
 - μ CT45:
 - 7.4 μ m, 1 stack, 1031 slices (7.64mm); Scan time: 48 mins
 - Analysis region
 - 100-200 slices, 0.5-1mm from the growth plate
 - Outcome measures: BV/TV, Tb.Th, Tb.N, Tb.Sp, SMI, Conn.D, BMD, TMD





Phenotyping

- Mouse cortical bone (long bone)
 - Midshaft of tibia or femur
 - μ CT35:
 - 6-10 μ m, 1 stack, ~210 slices (1.4-2.3mm); Scan time: 24 mins
 - μ CT45:
 - 7.4-10.4 μ m, 1 stack, 1031 slices (7.6-10.7mm); Scan time: 48 mins
 - Analysis region
 - Middle 50 slices
 - Outcome measures: Ct.Area, Ct.Th, pMOI, Ct.Po, TMD





Phenotyping

- Rat trabecular bone (long bone)
 - Proximal tibia
 - μ CT35:
 - 6-10 μ m, 2 stacks, ~420 slices (2.8-4.6mm); Scan time: 48 mins
 - μ CT45:
 - 7.4-10.4 μ m, 1 stack, 1031 slices (7.6-10.7mm); Scan time: 48 mins
 - Analysis region
 - ~200 slices, 1-2mm distal to the growth plate
 - Outcome measures: BV/TV, Tb.Th, Tb.N, Tb.Sp, SMI, Conn.D, TMD



Phenotyping

- Rat cortical bone (long bone)
 - Midshaft of tibia or femur
 - μ CT35:
 - 6-10 μ m, 1 stack, ~210 slices (1.4-2.3mm); Scan time: 24 mins
 - μ CT45:
 - 7.4-10.4 μ m, 1 stack, 1031 slices (7.6-10.7mm); Scan time: 48 mins
 - Analysis region
 - Middle 50 slices
 - Outcome measures: Ct.Area, Ct.Th, pMOI, Ct.Po, TMD



Phenotyping

- Mouse/Rat vertebral trabecular bone

- Lumbar vertebra L1-L4

- μ CT35:

- 6 μ m (mouse) / 6-10 μ m (rat)
- 2-3 stack, 420-630 slices
- Scan time: 48-72 mins

- μ CT45:



- 7.4-10.4 μ m, 1 stack, 1031 slices
- Scan time: 48 mins

- Scan Region: Between end plates

- Analysis region




- Between two end plates,
- Middle 150-200 slices

- Outcome measures: BV/TV, Tb.Th, Tb.N, Tb.Sp, SMI, Conn.D, BMD, TMD




Fracture Healing

- μ CT35:
 - 6-10 μ m
 - Scan vertically: 1-2 hours
 - Scan horizontally: 15-30 mins
- μ CT45:
 - 7.4-10.4 μ m
 - Scan vertically: 48 mins
 - Scan horizontally: 48 mins
- Scan region: Whole bone
- Analysis region: whole healing region excluding cortical and trabecular bone
- Outcome measures
 - BV: Callus Volume
 - Avg callus area = Callus Volume/ Total Length






Scaffolds

- PLA, PGA, PCL, etc.
 - Scan Region: Whole scaffold
 - μ CT35:
 - 6-10 μ m
 - Scan time: 0.5-2 hours
 - μ CT45:
 - 7.4-10.4 μ m
 - Scan time: 48 mins
 - Analysis Region: Whole scaffold
 - Outcome measure:
 - Total mineralized tissue content ($BV * TMD$)
 - Total mineralized tissue volume (BV)



Courtesy of Dr. Masahiro and Dr. Kenta Uchibe

μ CT 45 with a 20-sample holder carousel



Please write down the position# once you load the sample tube!

Note: Sample tubes from μ CT 35 are NOT compatible with μ CT 45.

μ CT 45 with a 20-sample holder carousel



Use the flashlight (we provided) to view sample placement

μ CT 45 with a 20-sample holder carousel

Demo of sample loading on uCT 45 <https://www.youtube.com/watch?v=NaxH4ycrIDg>



μ CT 45 with a 20-sample holder carousel

SCANCO MEDICAL
 μ CT 45

X-Ray
55 kVp
145 μ A

Carousel Position: **Click here after closing the door**

1 2 3 4 5 6 7 8 9 10 11 12 13 14 15 16 17 18 19 20

Check/Refresh **Unload** sample 5 is ready for measurement

Measurement Sample:

Number: **Click here if you want to open the door**

Name: Test Scan

Date of Birth:

Meas. Remarks:

Other... **Modify...** **New...**

Measurement History: *= aborted

1: No.:	2	CTR:	4	100 Slices	*	11-NOV-2019	10:37 P
2: No.:	29	CTR:	12	1032 Slices		2-DEC-2019	16:26
3: No.:	30			633 Slices		2-DEC-2019	17:44

Controlfile:

4:	QC1 90kVp, 0.1mm Cu, 1200mg HA (custom)
5:	QC2 70kVp, 0.5mm Al, 1200mg HA (monthly)
10:	Test Scan
11:	PCMD-55kVP-3 μ m-300ms-FOV10.2mm
12:	PCMD-55kVP-4.5 μ m-300ms-FOV15.3mm
14:	PCMD-55kVP-7.4 μ m-300ms-FOV25.2mm
15:	PCMD-55kVP-10.4 μ m-300ms-FOV35.4mm
16:	PCMD-55kVP-14.6 μ m-300ms-FOV49.6mm


Modify... **Delete...** **New...** **Lock...**

Scout-View... **Scan...** **Reset**

Add Scout **Add Scan** **Task List...**



Please confirm

 Insert or remove the sample holders now!
To start a measurement, please close the door and press yes.

Do you want to start a measurement?

Yes **No** **Cancel**

Scout-View for multiple tubes

(1)

1 2 3 4 5 6 7 8 9 10 11 12 13 14 15 16 17 18 19 20

Check/Refresh Unload sample 5 is ready for measurement

Measurement Sample:

Number: 3

Name: Test Scan

Date of Birth:

Meas. Remarks:

Other... Modify... New...

Measurement History: * = aborted

1: No.:	2	CTR:	4	100 Slices	* 11-NOV-2019 10:37 P
2: No.:	29	CTR:	12	1032 Slices	2-DEC-2019 16:26
3: No.:	30			633 Slices	2-DEC-2019 17:44

Controlfile:

```

4: QC1 90kVp, 0.1mm Cu, 1200mg HA (custom)
5: QC2 70kVp, 0.5mm Al, 1200mg HA (monthly)
10: Test Scan
11: PCMD-55kVP-3um-300ms-FOV10.2mm
12: PCMD-55kVP-4.5um-300ms-FOV15.3mm
14: PCMD-55kVP-7.4um-300ms-FOV25.2mm
15: PCMD-55kVP-10.4um-300ms-FOV35.4mm
(2) 16: PCMD-55kVP-14.6um-300ms-FOV49.6mm
    
```

Modify... Delete... New... Lock...

Scout-View... Scan... Reset

(3) Add Scout Add Scan Task List... (4)

50 0.000

Tasks

- Done	Scout-View S:	3, P: 5 C:	12 #:1032 H: 9	3-DEC-2019 09:40:37
- Done	Scout-View S:	3, P: 5 C:	14 #:1031 H: 8	3-DEC-2019 09:40:38
- Done	Scout-View S:	3, P: 5 C:	15 #:1031 H: 7	3-DEC-2019 09:40:39
- Done	Scout-View S:	3, P: 5 C:	16 #:1031 H: 6	3-DEC-2019 09:40:41
1: Pending	CT-Scan S:	3, P: 5 C:	12 #:1032 H: 9	3-DEC-2019 10:05:31
- Pending	CT-Scan S:	3, P: 5 C:	14 #:1031 H: 8	3-DEC-2019 10:05:55
- Pending	CT-Scan S:	3, P: 5 C:	15 #:1031 H: 7	3-DEC-2019 10:06:14
- Pending	CT-Scan S:	3, P: 5 C:	16 #:1031 H: 6	3-DEC-2019 10:06:29

Remove... Active Keep Load Scout (6)

Current Task:

3: Test Scan
 CT-Scan: ø 49.640 mm, 61.281-76.334 mm, H: U45824, ø 48mm x H 110mm
 55 kVp, 145 µA, AL 0.5mm, 1x300 ms
 Dim: 3400x3400, 1031 slices, Voxelsize 14.6 µm

(5) Start Interact. Tasks Submit Batch Scans Clear List Close Window

Scout-View for multiple tubes

- We suggest you do scout view for all tubes first before selecting scanning regions:

1. Select carousel position#

2. Select Controlfile#

3. Click "Add Scout"

Do NOT click "Scout-View..."

4. All scout tasks are listed in the "Task List..."

Carousel Position:

1 2 3 4 5 6 7 8 9 10 11

Check/Refresh Unload found 3 valid samples on carousel

Controlfile: []

4: QC1 90kVp, 0.1mm Cu, 1200mg HA (custom)
5: QC2 70kVp, 0.5mm Al, 1200mg HA (monthly)
10: Test Scan
11: PCMD-55kVP-3um-300ms-FOV10.2mm
12: PCMD-55kVP-4.5um-300ms-FOV15.3mm
14: PCMD-55kVP-7.4um-300ms-FOV25.2mm
15: PCMD-55kVP-10.4um-300ms-FOV35.4mm
16: PCMD-55kVP-14.6um-300ms-FOV49.6mm

Modify... Delete... New... Lock...

Scout-View... Scan... Reset

Add Scout Add Scout Task List...


Tasks

- Done	Scout-View S:	3, P: 5 C:	12 #:1032 H: 9	3-DEC-2019 09:40:37
- Done	Scout-View S:	3, P: 5 C:	14 #:1031 H: 8	3-DEC-2019 09:40:38
- Done	Scout-View S:	3, P: 5 C:	15 #:1031 H: 7	3-DEC-2019 09:40:39
- Done	Scout-View S:	3, P: 5 C:	16 #:1031 H: 6	3-DEC-2019 09:40:41
1: Pending	CT-Scan S:	3, P: 5 C:	12 #:1032 H: 9	3-DEC-2019 10:05:31
- Pending	CT-Scan S:	3, P: 5 C:	14 #:1031 H: 8	3-DEC-2019 10:05:55
- Pending	CT-Scan S:	3, P: 5 C:	15 #:1031 H: 7	3-DEC-2019 10:06:14
- Pending	CT-Scan S:	3, P: 5 C:	16 #:1031 H: 6	3-DEC-2019 10:06:29

Remove... Active Keep Load Scout

Scout-View for multiple tubes


- We suggest you do scout view for all tubes first before selecting scanning region:
 5. Repeat steps 1~3 until all tubes are in “Task List...”, then click “Start Interact. Tasks”



Take a break. Scout view may take ~5 min/tube.

(If you have 10 tubes, you can come back in 1 hour)

5. After scout view for all tubes are done, load each scout view (“Load Scout”) to select scanning region.



Scout-View for multiple tubes

Tasks

-

Done

Scout-View

S:

3,

P:

5

C:

12

#:

1032

H:

9

3-DEC-2019

09:40:37

-

Done

Scout-View

S:

3,

P:

5

C:

14

#:

1031

H:

8

3-DEC-2019

09:40:38

-

Done

Scout-View

S:

3,

P:

5

C:

15

#:

1031

H:

7

3-DEC-2019

09:40:39

-

Done

Scout-View

S:

3,

P:

5

C:

16

#:

1031

H:

6

3-DEC-2019

09:40:41

1:

Pending

CT-Scan

S:

3,

P:

5

C:

12

#:

1032

H:

9

3-DEC-2019

10:05:31

-

Pending

CT-Scan

S:

3,

P:

5

C:

14

#:

1031

H:

8

3-DEC-2019

10:05:55

-

Pending

CT-Scan

S:

3,

P:

5

C:

15

#:

1031

H:

7

3-DEC-2019

10:06:14

-

Pending

CT-Scan

S:

3,

P:

5

C:

16

#:

1031

H:

6

3-DEC-2019

10:06:29

▲

▼

Remove...

☐ Active

☐ Keep

Load Scout

Current Task:

3: Test Scan

CT-Scan: ø 49.640 mm, 61.281–76.334 mm, H: U45824, ø 48mm x H 110mm

55 kVp, 145 µA, AL 0.5mm, 1x300 ms

Dim: 3400x3400, 1031 slices, Voxelsize 14.6 µm

Start Interact. Tasks

Submit Batch Scans

Clear List

Close Window

The logo of the University of Pennsylvania, featuring a shield with a red top section containing three open books, a blue bottom section containing three white stars, and a white diagonal section containing a red anchor.

Penn

Scout-View for multiple tubes

Tasks

		Sample#	Controlfile#							
- Done	Scout-View	S: 3	P: 5	C: 12	#:1032	H: 9	3-DEC-2019	09:40:37		
- Done	Scout-View	S: 3	P: 5	C: 14	#:1031	H: 8	3-DEC-2019	09:40:38		
- Done	Scout-View	S: 3	P: 5	C: 15	#:1031	H: 7	3-DEC-2019	09:40:39		
- Done	Scout-View	S: 3	P: 5	C: 16	#:1031	H: 6	3-DEC-2019	09:40:41		
1: Pending	CT-Scan	S: 3	P: 5	C: 12	#:1032	H: 9	3-DEC-2019	10:05:31		
- Pending	CT-Scan	S: 3	P: 5	C: 14	#:1031	H: 8	3-DEC-2019	10:05:55		
- Pending	CT-Scan	S: 3	P: 5	C: 15	#:1031	H: 7	3-DEC-2019	10:06:14		
- Pending	CT-Scan	S: 3	P: 5	C: 16	#:1031	H: 6	3-DEC-2019	10:06:29		

Position#

Slice#

▲

▼

Remove...

☐ Active

☐ Keep

Load Scout

Current Task:

3: Test Scan
CT-Scan: ø 49.640 mm, 61.281–76.334 mm, H: U45824, ø 48mm x H 110mm
55 kVp, 145 µA, AL 0.5mm, 1x300 ms
Dim: 3400x3400, 1031 slices, Voxelsize 14.6 µm

Start Interact. Tasks

Submit Batch Scans

Clear List


Close Window

Room 335A: μ CT35 & μ CT45



Penn

Room 335A: μ CT45



Penn

Recent Update: Basic tools are provided



- Basic tools, ethanol, PBS, and markers are provided
- Allowing users to wear gloves on waterproof keyboard
- Please wipe keyboard with 70% ethanol after use



Penn

Room 368C: vivaCT 40



Recent Update: Basic tools are provided



- Basic tools, ethanol, PBS, and markers are provided
- Allowing users to wear gloves on waterproof keyboard
- Please wipe keyboard with 70% ethanol after use




Penn

Room 368B: vivaCT 80




Room 368B: vivaCT 80



Penn

Room 315: Analysis Computer



- Before use, restart the computer to avoid previous errors
- For calendar access, send request to pcmd.microct@gmail.com




Penn

Room 315: Analysis Computer

- Video tutorial <https://www.youtube.com/watch?v=Fz0KfaLYDKw>
(You may adjust your preferred video speed in Youtube)



Room 315: Analysis Computer



- Policy: No hours limit, but 30-min “No Show” is NOT allowed.
Calendars for analysis computers:
[MicroCT Analysis \(I\) - Room 315](#)
[MicroCT Analysis \(II\) - Room 315](#)



Penn

Booking Policies for μ CT Scanner

- Calendars for scanners:

[MicroCT35](#)

[MicroCT45](#)

[VivaCT40](#)

[VivaCT80](#)

- Label your Google event “PI’s initial_PennKey_scan only_[Phone Number](#)”, e.g.:
XSL_yiluz_scan & analysis_1234567890
- 30 min “No Scan” is NOT allowed. Adjust your reservation if your plans change!**
- No restrictions for off-hours (< 9am or > 5pm) & weekend.
- No restrictions for *in vivo* (live animal) scanning.
- Currently, no restrictions for vivaCT 80 and μ CT 45.
- Scanning user always has priority over “analysis only” user.
- μ CT 35 & vivaCT 40:
4 hours maximum per lab per weekday ([9-1pm or 1-5pm](#))
(Exception: On your scanning day, you can reserve if additional time is available)



Scanning Policies

Sample naming: **Starts with PI's initial** + user's initial or first name or pennkey (for billing purposes)

Scanning:

- All *in vivo* μ CT scans require IACUC approval
- In the case of error messages, please capture photos of the error message, and send to pcmd.microct@gmail.com
- Any scan extending **10 min** into other user's period will be stopped by manager

After CT scanner use:

- Clean up
- Do NOT turn off the CT scanner
- Do NOT “session out” the system
- **Samples will be trashed if left over 2 days after scanning**

Data Storage Policies

- Google Drive auto-share: Requested files will be removed (from Google drive) 14 days after sharing
- Image files will be removed from servers to offline tapes within 3 months after scanning.
(If you are still analyzing samples scanned 3 months ago, please let us know in advance!)
- Retrieving image files from tapes is very time consuming and depends on the availability of both tape drive and server space.

Please request tape file retrieval at least 7 days in advance!


microCT Sample Database

MicroCT_Sample_Database

File Edit View Insert Format Data Tools Add-ons Help

100% View only

	A	B	C	D	E	F	G	H	I	J	K
1	2019-11-26 16:03 Last updated by MicroCT35		MicroCT35	sorted by scan date		File & Tape Request Form (v201912)					
2	PI User	Sample Name	Scan starts	Scan ends	Sample#	Measure#	RSQ Size	ISQ Size	File Types	Analyzed?	Controlfile
3	FXL Fangfang	Long_Fang_8511	2019-11-26 18:07	18:32	7527	18383					109: PCMD-bone-6um-D11.5mm-300ms
4	FXL Fangfang	Long_Fang_8511	2019-11-26 17:19	18:07	7527	18382					109: PCMD-bone-6um-D11.5mm-300ms
5	FXL Fangfang	Long_Fang_8533	2019-11-26 16:31	17:19	7526	18381					109: PCMD-bone-6um-D11.5mm-300ms
6	FXL Fangfang	Long_Fang_8533	2019-11-26 16:06	16:31	7526	18380					109: PCMD-bone-6um-D11.5mm-300ms
7	FXL Fangfang	Long_Fang_8509	2019-11-26 15:41	16:06	7525	18379	0KB				109: PCMD-bone-6um-D11.5mm-300ms
8	FXL Fangfang	Long_Fang_8509	2019-11-26 14:52	15:40	7525	18378	1.84GB	3.15GB			109: PCMD-bone-6um-D11.5mm-300ms
9	FXL Fangfang	Long_Fang_8534	2019-11-26 14:04	14:52	7524	18377	1.84GB	3.50GB			109: PCMD-bone-6um-D11.5mm-300ms
10	FXL Fangfang	Long_Fang_8534	2019-11-26 13:31	13:56	7524	18376	969.72MB	1.79GB			109: PCMD-bone-6um-D11.5mm-300ms
11	SY Yang	sy-ly-r11m20tazoc-3m	2019-11-26 8:51	9:39	7523	18375	1.81GB	3.46GB			10: Test Scans
12	SY Yang	sy-ly-r5m20tazoc-3m	2019-11-26 8:03	8:51	7522	18374	1.81GB	3.46GB			10: Test Scans
13	SY Yang	sy-ly-367r1moc-3m	2019-11-26 7:09	7:57	7521	18373	1.81GB	3.46GB			10: Test Scans
14	FXL Chao	Long lab_chao_7483 N2P	2019-11-26 3:39	5:10	7520	18372	3.73GB	7.09GB			109: PCMD-bone-6um-D11.5mm-300ms
15	FXL Chao	Long lab_chao_7482 N2p	2019-11-26 2:08	3:39	7519	18371	3.73GB	7.09GB			109: PCMD-bone-6um-D11.5mm-300ms
16	FXL Chao	Long lab_chao_7477 N2P	2019-11-26 0:37	2:08	7518	18370	3.73GB	7.09GB			109: PCMD-bone-6um-D11.5mm-300ms
17	FXL Chao	Long lab_chao_7467 N2P	2019-11-25 22:59	11-26 0:30	7517	18369	3.73GB	7.09GB			109: PCMD-bone-6um-D11.5mm-300ms
18	FXL Chao	Long lab_chao_7599 N2P	2019-11-25 21:11	22:42	7516	18368	3.73GB	7.09GB			109: PCMD-bone-6um-D11.5mm-300ms
19	FXL Chao	Long lab_chao_7476 N2P	2019-11-25 19:40	21:11	7515	18367	3.73GB	7.09GB	AIM,GOBJ,TXT	Analyzed	109: PCMD-bone-6um-D11.5mm-300ms
20	FXL Chao	Long lab_chao_7460 N2P	2019-11-25 18:08	19:39	7514	18366	3.73GB	7.09GB			109: PCMD-bone-6um-D11.5mm-300ms
21	FXL Chao	Long lab_chao_7451 N2P	2019-11-25 16:37	18:08	7513	18365	3.73GB	7.09GB			109: PCMD-bone-6um-D11.5mm-300ms


CT35_By_date
Viva40_By_date
Viva80_By_date
CT45_By_date
CT35_By_sample
Viva40_By_sample
Viva80_By_sample

Sort by scan date

Sort by sample#

Automated File Request

MicroCT_Sample_Database

File Edit View Insert Format Data Tools Add-ons Help

100% View only

	A	B	C	D	E	F	G	H	I	J	K
1	2019-11-26 16:03 Last updated by MicroCT35		MicroCT35	sorted by scan date			File & Tape Request Form (v201912)				
2	PI User	Sample Name	Scan starts	Scan ends	Sample#	Measure#	RSQ Size	ISO Size	File Types	Analyzed?	Controlfile
3	FXL Fangfang	Long_Fang_8511	2019-11-26 18:07	18:32	7527	18383					109: PCMD-bone-6um-D11.5mm-300ms

- Fill in this file request form, and send to pcmd.microct@gmail.com

Your Gmail: meniscus@gmail.com								
MicroCT35			Vivact40			Vivact80		
Sample#	Measure#	File_Types	Sample#	Measure#	File_Types	Sample#	Measure#	File_Types
7234	17817	DICOM						
7234	17816	DICOM						
7234	17815	DICOM						
7234	17814	DICOM						
7234	17813	DICOM						
7234	17812	DICOM						
7233	17811	DICOM						
7233	17810	DICOM						
7233	17809	DICOM						
7204	17726	DICOM						
7204	17725	DICOM						
7204	17724	DICOM						

Example (1)



Penn

Automated File Request

- Fill in this file request form, and send to pcmd.microct@gmail.com

Your Gmail: meniscus@gmail.com , tendon@gmail.com , muscle@gmail.com , bone@gmail.com								
MicroCT35			Vivact40			Vivact80		
Sample#	Measure#	File_Types	Sample#	Measure#	File_Types	Sample#	Measure#	File_Types
						134	839	AIM, TXT
						134	840	AIM, TXT
						134	841	AIM, TXT
						134	842	AIM, TXT
						134	843	AIM, TXT
						134	844	AIM, TXT
						134	845	AIM, TXT
						134	846	AIM, TXT
						134	847	AIM, TXT
						135	850	AIM, TXT
						136	851	AIM, TXT

Example (2)

After analysis, we suggest you keep a copy of GOBJ and TXT.
Analysis results are in the TXT file with the name "3DRESULTS".
(Open it with Excel, NOT notepad.)

Tape Retrieval Request

- If your samples have been moved into tapes, please select **"YES (They are in tape)"** at cell P2
- Send to pcmd.microct@gmail.com
- Tape retrieval may take over 1 week!
(Please send to us as early as possible)

Are you requesting from tapes?	
NO (default)	<input type="radio"/>
YES (They are in tape)	<input checked="" type="radio"/>

please select "YES (They are in tape)" at cell P2. (See [Example 3](#))



Vivact80			MicroCT45		
Sample#	Measure#		Sample#	Measure#	
260	940				
261	939				
262	944				
263	943				
264	942				
265	947				
266	946				

**Example (3)
For Tape Retrieval**

Are you requesting from tapes?
YES (They are in tape)

Automated System for Sample Re-alignment

- Examples:
 - Scanning angle is too far away
 - From transverse to sagittal or coronal plane



1	2. Fill in all yellow and red boxes (4 boxes: B2, A4, B4, C4)					
2	3. Fill in this Excel sheet, 1 row for ONLY 1 measure number					
3	4. Send this Excel sheet to pcmd.microct@gmail.com					
4	1. Enter your email here: bone@gmail.com					
5	2. Select machine:	3. Select rotating options:	4. Select import options:			
6	MicroCT35	3 Point Rotation (AlignZ)	Import with new sample name	Use COMMA btw each number, e.g.: 100,100,100		
7	Sample#	Measure#	New Sample#	P1	P2	P3
8	1606	8335	1608	174,224,0	173,99,0	93,154,0
9						
10						

Current status: In development (expected 02/2020)



Penn

Questions?