

## *PCMD MicroCT Imaging Core Learning Lunch Series*

# Deep Learning for MicroCT Image Segmentation via Dragonfly

April 28<sup>th</sup>, 2022

PCMD MicroCT Imaging Core

(Management team: X. Sherry Liu, Yilu Zhou)

## Invited speakers:

Dr. Ryan Locke from Dr. Robert Mauck's group (Penn McKay Lab)

Ranhui Xi from Dr. Marco Tizzano's group (Penn Dental)



PENN

CENTER for

MUSCULOSKELETAL

DISORDERS

Penn Center for Musculoskeletal Disorders (P30-AR069619)



# Outlines

- Brief introduction of our core facility
- Recent progress of our video tutorials, automated services
- The Dragonfly software & our new workstation
- Dragonfly application highlights from Dr. Ryan Locke  
(Dr. Robert Mauck's group, Penn McKay Lab)
- Dragonfly application highlights from Ranhui Xi  
(Dr. Marco Tizzano's group, Penn Dental)
- Q & A



# μCT Imaging Core Resources

	Model	Location	Scan Size (ØxL;mm)	Voxel Size (μm)	Applications
1	μCT 35	Stemmler Hall	37.9 x 120	3.5-72	High resolution <i>ex vivo</i> scans
2	μCT 45	Stemmler Hall	50 x 120	3.0-100	High resolution <i>ex vivo</i> scans
3	vivaCT 40	Stemmler Hall	38.9 x 145	10.5-76	High resolution <i>in vivo</i> scans for small animals
4	vivaCT 80	Stemmler Hall	80 x 145	10.4-76	High resolution <i>in vivo</i> scans for small animals
5	μCT 50	PVAMC/TMRC	50 x 120	0.5-100	Ultra high resolution (sub-micron) <i>ex vivo</i> scans
6	vivaCT 75	PVAMC/TMRC	79.9 x 145	21-150	<i>In vivo</i> scans for small animals; <i>ex vivo</i> scans for large specimens
7	XtremeCT II	CTRC	140 x 200	60-82	Clinical scans for peripheral skeleton


# Ex vivo (Specimen) Scanners

- $\mu$ CT 35 (Purchased in 2012)
  - Voxel sizes: 3.5  $\mu$ m, 6  $\mu$ m, 10  $\mu$ m, 15  $\mu$ m, 18.5  $\mu$ m



# Ex vivo (Specimen) Scanners

- $\mu$ CT 45 (Purchased in 2019 *new!*)
  - Voxel sizes: 3  $\mu$ m, 4.5  $\mu$ m, 7.4  $\mu$ m, 10.4  $\mu$ m, 14.6  $\mu$ m
  - Carousel system supporting 20 sample holders
  - “Air” filter for scanning low density materials
  - “Copper” filter for scanning metal implant





# In vivo Scanners

- vivaCT 40 (Purchased in 2010)
  - Voxel sizes: 10.5  $\mu\text{m}$ , 12.5  $\mu\text{m}$ , 15  $\mu\text{m}$ , 17.5  $\mu\text{m}$ , 19  $\mu\text{m}$




# In vivo Scanners

- vivaCT 80 (Purchased in 2018 *new!*)
  - Voxel sizes: 10.4  $\mu\text{m}$ , 11.6  $\mu\text{m}$ , 13  $\mu\text{m}$ , 16.1  $\mu\text{m}$ , 20.8  $\mu\text{m}$ , 26  $\mu\text{m}$
  - Internal heating device to keep animal warm
  - Camera to monitor animal's breathing
  - Ex vivo scan for specimen from large animals or human cadaver



# MicroCT Analysis PC

- 2 PCs for MicroCT Analysis (315 Stemmler)
  - Windows 10 platform
  - Either remote or onsite access
  - Scanco software





# Video Tutorials & Instruction Documents

<https://www.med.upenn.edu/pcmd/mctimagingcore/user-tutorials.html>

<https://www.youtube.com/channel/UCzznR9Fdv-3kjEX7miwsi0A>

## Video Tutorials:

Our YouTube channel: <https://www.youtube.com/channel/UCzznR9Fdv-3kjEX7miwsi0A/>

### μCT scan setup:

1. How to set up a scan on μCT35 ([PDF download](#)) ([Video download](#))

<https://www.youtube.com/watch?v=QUtoQqIYJ8o>

2. Demo: How to set up a scan on μCT45 (Recommended: Carousel version)

([PDF download](#)) ([Video download](#))

Note: To use this Carousel version, please remove the sample holder on the rotation stage.

<https://www.youtube.com/watch?v=fz1ffR5XyE>

3. Demo: How to set up a scan on μCT45 (Non-carousel version) ([PDF download](#)) ([Video download](#))

Note: To use this Non-carousel version, please remove all sample holders on the carousel.

<https://www.youtube.com/watch?v=JEoLn1igEjE>

4. How to set up an ex vivo scan on VivaCT40 ([PDF download](#)) ([Video download](#))

<https://www.youtube.com/watch?v=sxvTV4bvosw>

5. How to set up an ex vivo scan on VivaCT80 ([PDF download](#)) ([Video download](#))

<https://www.youtube.com/watch?v=HdQYVxwjuIXM>

### μCT viewing & analysis:

1. How to use "microCT Analysis" computers ([PDF download](#)) ([Video download](#))

<https://www.youtube.com/watch?v=qHHcB6KJJe4>

2. Tutorial for cropping, exporting, and requesting microCT images ([PDF download](#)) ([Video download](#))


<https://www.youtube.com/watch?v=umRF6ODcQqQ>

3. Tutorial for 3D display of microCT images ([PDF download](#)) ([Video download](#))

<https://www.youtube.com/watch?v=YdQSO41rgR8>

4. Tutorial for cortical bone analysis (mouse tibia midshaft) ([PDF download](#)) ([Video download](#))

<https://www.youtube.com/watch?v=B4OEgX8Bkwg>



culoskeletal Disorders



# Publications from our users

- Our users have published over 230 journal articles on their  $\mu$ CT projects.
- Selected publications with detailed  $\mu$ CT protocols for other users to cross reference:

<https://www.med.upenn.edu/pcmd/mctimagingcore/publications.html>

## 1. Calcified Tissue Imaging

### 1.1. Skeletal Phenotyping

#### 1.1.1. Rodents

OA study (gene therapy): Proximal tibia of Sprague-Dawley rats. Mason, J.B., et al., *Wnt10b and Dkk-1 gene therapy differentially affects bone mass and osteophytosis in a skeletally mature rat model of osteoarthritis*. J Orthop Res, 2016. 34(10): p. 1947-1958.

OA study (DMM Model): Medial epiphysis of the mice. Sambamurthy, N., et al., *Chemokine receptor-7 (CCR7) deficiency exacerbates bone loss and osteoarthritis in a murine model of osteoarthritis*. J Orthop Res, 2016. 34(10): p. 1959-1970.

## 2. Non-calcified Tissue Imaging

### 2.1. Cartilage Imaging

Cartilage repair: Osteochondral specimens from mice. Friedman, J.M., et al., *Comparison of Fixation Techniques for Small Animal Cartilage Repair Models*. Cartilage, 2017. 28(3): p. 375-385.

Weightbearing Porcine Large Animal Model. Cartilage repair. Pfeifer, C.G., et al., *Age-Dependent Subchondral Bone Loss in a Porcine Model of Osteoarthritis*. Part C Methods, 2017. 23(11): p. 745-753.

Patel, J.M., et al., *Resorbable pins to enhance scaffold integration in a murine model of osteoarthritis*. PLoS One, 2017. 12(3): p. e0173995.

## 3. In Vivo Small Animal Imaging

Reproducibility and Radiation study: Mice distal femur. Zhao, H., et al., *Reproducibility and Radiation Effect of  $\mu$ CT Imaging in Mice*. Ann Biomed Eng, 2016. 44(10): p. 1000-1010.

Bone remodeling study: Longitudinal *in vivo* scanning. de Bakker, C.M.J., et al., *Minimizing Interpolation Bias in Longitudinal  $\mu$ CT Scans*. Ann Biomed Eng, 2016. 44(10): p. 1011-1020.

Reproduction cycles study: Longitudinal *in vivo* scanning. de Bakker, C.M., et al., *Adaptations in the Microarchitectural Response to Multiple Reproductive Cycles in Rats*. J Bone Miner Res, 2017. 32(1): p. 1-12.

## 4. Clinical Imaging

HR-pQCT scanner (XtremeCT II), human tibia. Zhao, X., et al., *Feasibility of assessing bone material properties using HR-pQCT*. PLoS One, 2017. 12(3): p. e0173995.

## 5. Other Imaging

Metal implants in rat brain (90 kVp with a copper filter). Burton A, et al., *Wireless, battery-free, and fully implantable  $\mu$ CT for in vivo imaging*. Nanoeng. 2021;7:62.



# Automated Services

- File request: fully automated service sharing MicroCT files to users (running 7/24)
- Auto compiling microCT results into Excel sheet

[https://www.med.upenn.edu/orl/uct/assets/user-content/secure/User\\_file\\_request%20\(v2020.01\).xlsx](https://www.med.upenn.edu/orl/uct/assets/user-content/secure/User_file_request%20(v2020.01).xlsx)

- Sample realignment/reorientation request: fully automated service help users to do sample realignment (running 7/24)

[https://www.med.upenn.edu/orl/uct/assets/user-content/secure/Sample\\_Realignment\\_request\(v2020.01\).xlsx](https://www.med.upenn.edu/orl/uct/assets/user-content/secure/Sample_Realignment_request(v2020.01).xlsx)



# Dragonfly Software

- Free non-Commercial License <https://www.theobjects.com/dragonfly/get-non-commercial-licensing-program.html>
- Supports many image formats: dcm, tif, Scanco aim, czi, raw, pic, img, etc...
- Training video tutorials available <https://www.theobjects.com/dragonfly/tutorials.html>

DRAGONFLY  
welcome


FEATURES  
key benefits

DEEP LEARNING  
advanced processing

RESOURCES  
learn more

SUPPORT  
get help

GET DRAGONFLY  
achieve real results

Search... 

Resources / Videos and Tutorials / Dragonfly Tutorials

Understanding the Dragonfly interface and its wide range of features is essential to improving your workflow.

Additional videos featuring in-depth lessons on Dragonfly are available in the section [In-Depth Lessons](#).

All

Key Features


Image Processing

Segmentation/Quantification

Presentation


Specialized Applications

Macros




Getting Started

Keywords:  
UI, panel, scene, view, 2D, 3D, state, action, data properties, workspace



Importing Images

Keywords:  
TIFF, RAW, image settings, volumetric, stack, crop, transform, offset, slope



Importing DICOM Images

Keywords:  
DICOM, attributes, study, series, preferences, server, folder



# Dragonfly Workstation

- New workstation for Dragonfly software (324 Stemmler)
  - Windows 10 platform
  - PMACS account required (either remote or onsite access)
  - Hard drive: 18 TB
  - Memory: 128 GB
  - GPU: NVIDIA RTX 3090 (Ready in June)

