



## Hybridoma cell culture and antibody purification

1. Thawing Hybridoma cells
  - a. Cells are more delicate than normal tissue culture cells
  - b. Thaw cells rapidly, then transfer them to 5 ml warm RPMI 1640 medium containing 10% FBS and 1:100 P/S (Mediatech MT30002CI)
  - c. Spin down at 800 x g for 5 min
  - d. Remove supe and resuspend gently in 5 ml warm RPMI 1640
  - e. Plate onto 6cm tissue culture plate
  - f. If low on vials of cells, freeze some down in freezing medium
  
2. Splitting cells
  - a. Cells need to be split 1:5 every 3 days or so.
    - i. This is usually done when medium turns from pink to orange
      1. Density should be around  $1 \times 10^6$  cells/ml when splitting
      2. Do not let cells get overgrown
    - ii. Cells are not adherent, should look round, smooth and quite large
  - b. Split cells into larger plates every-time by simply adding cells to fresh medium
    - i. 5 ml into 20 on large 15 cm plates
    - ii. 25 ml into 100 ml in 75 cm<sup>2</sup> flasks
    - iii. 125 ml into 5x 100 ml in 5x 75 ml cm<sup>2</sup> flasks
  - c. Never throw away any supernatant as it contains the secreted antibody
  
3. Antibody Purification
  - a. Collect all cells, ~500 ml, and spin down at 2.5k, 5 min 4 °C
    - i. Collect supernatant and filter with 0.22 µm filter top
  - b. Equilibrate 1 ml HiTrap Protein G HP column (GE healthcare 17-0404-03) in cold room
    - i. Cold, filtered 10 ml PBS
  - c. Apply supernatant to column at 1.5 ml/min
  - d. Wash column with 20 ml PBS
  - e. Elute antibodies with 5 ml 100 mM Glycine pH 2.5
    - i. Collect 0.5 ml fractions into tubes containing 25 µL 1M TRIS pH 9.5
    - ii. Gently mix tubes immediately after elution is complete to neutralize pH
  - f. Wash column with 20 ml PBS
    - i. Store in 20% EtOH at 4 °C
  - g. Run fractions on SDS PAGE gel
    - i. Should get two bands, 30 & 50 kDa
    - ii. Some samples may be yellow due to pH, they will turn blue when loaded onto gel

- iii. Can also check by abs@ 280nm, but since only 10 samples, you will run them all on a gel anyway
- h. Collect fractions containing high concentrations of antibody
  - i. Dialyze in 2 L PBS over night
  - ii. If time permits, change buffer once to remove excess glycerol
- i. Dialyze in 500 ml of PBS/50% Glycerol for at least 5 hours
- j. Measure antibody concentration
  - i. Should get ~300  $\mu$ Ls or 5.0 mg/ml
- k. Store at -20 °C