


NADH Enzyme-Linked Assay



Stock Solutions

4000 U/ml	Lactic Dehydrogenase (Sigma L-1254), store in freezer in 50% glycerol
10,000 U/ml	Pyruvate Kinase (Sigma P-9136), store in freezer in 50% glycerol
100 mM	Phospho(enol)pyruvate (Sigma, P-7127), store in freezer in water
1 mg	NADH, reduced form (Sigma, N0786), store at room temperature in the dark

5x Cocktail Solution

Add 1.2 ml assay buffer (usually KMg50 or KMg10) to 1 mg NADH, and add:

32 μ l	Lactic Dehydrogenase
64 μ l	Pyruvate Kinase
32 μ l	Phospho(enol)pyruvate

This will result in the following concentrations:

1000 μ M	NADH
100 U/ml	Lactic Dehydrogenase
500 U/ml	Pyruvate Kinase
2.5 mM	Phospho(enol)pyruvate

Assay

1. When using the stopped flow, one syringe should contain a 2x mixture of the Cocktail Solution and the appropriate 2x MgATP concentration (usually 4 mM).

Example of 2x Cocktail (2.88 ml):

1200 μ l	5x Cocktail from above
115 μ l	0.1 M ATP
11.5 μ l	1 M MgCl ₂
1550 μ l	Assay Buffer

2. Detect absorbance change at 340 nm, or detect fluorescence emission (Ex. 340 nm) with 400 nm long-pass filter.

3. Use extinction coefficient for NADH ($\epsilon_{340} = 6220 \text{ M}^{-1}\text{cm}^{-1}$) to determine rate of ADP production (NADH loss). Remember that the path length of the stopped-flow cuvette is 0.2 cm.

4. When using fluorescence detection, a standard curve with known ADP concentrations must be obtained.