

## ANTIBIOTICS

### Ampicillin

*Stock solution.* 25 mg/ml of the sodium salt of ampicillin in water. Sterilize by filtration and store in aliquots at  $-20^{\circ}\text{C}$ .

*Working concentration.* 35–50  $\mu\text{g/ml}$ .

### Chloramphenicol

*Stock solutions.* 34 mg/ml in 100% ethanol. Store at  $-20^{\circ}\text{C}$ .

*Working concentration.* For amplification of plasmids, 170  $\mu\text{g/ml}$ ; for selection of resistant bacteria, 30  $\mu\text{g/ml}$ .

### Kanamycin

*Stock solution.* 25 mg/ml in water. Sterilize by filtration and store in aliquots at  $-20^{\circ}\text{C}$ .

*Working concentration.* 50  $\mu\text{g/ml}$ .

### Nalidixic Acid

*Stock solution.* 20 mg/ml in water. Sterilize by filtration and store in aliquots at  $-20^{\circ}\text{C}$ .

*Working concentration.* 20  $\mu\text{g/ml}$ .

### Streptomycin

*Stock solution.* 20 mg/ml in water. Sterilize by filtration and store in aliquots at  $-20^{\circ}\text{C}$ .

*Working concentration.* 25  $\mu\text{g/ml}$ .

### Tetracycline

*Stock solution.* 12.5 mg/ml tetracycline hydrochloride in ethanol/water (50% v/v). Store at  $-20^{\circ}\text{C}$ .

*Note.* Because tetracycline is light-sensitive, solutions and plates containing the antibiotic should be stored in the dark.

*Working concentration.* 12.5–15.0  $\mu\text{g/ml}$ . Magnesium ions are antagonists of tetracycline. Use media without magnesium salts (e.g., LB) for selection of bacteria resistant to tetracycline.