

Draft

**Ca<sup>++</sup> Flux Analysis via FACS**

protocol by Susannah D. Barbee

1. Prepare single-cell suspensions of primary thymocytes or lymphocytes. Resuspend cells at 10<sup>6</sup>/mL in RT or 37°C Cell-Loading Buffer (CLB: HBSS + 1% FBS).
2. Add Indo-1 (2 mg/mL in DMSO, Molecular Probes I-1203) to a final concentration of 2 µg/mL.
3. Incubate 30 min 37°C.
4. For simultaneous staining of surface antigens: stain with azide-free FITC- and PE-conjugated antibodies for 20 min at RT<sup>+</sup>.
5. Resuspend in CLB at 5 x 10<sup>5</sup>/mL. Keep at RT until analysis (within 2 hrs).
6. Incubate sample aliquots for 5 min at 37°C before running.
7. Collect a baseline ratio for 30 sec before adding stimulants, collect 8-10 min.
8. Use Ionomycin to measure maximal flux: dilute stock (500 µg/mL in DMSO) 1/5 with CLB to 100 µg/mL. Add 10 or 20 µL per 1 mL sample for 1 or 2 µg/mL final concentration.  
\*\* Always flush lines with DMSO and then CLB after stimulations to remove any excess ionomycin.

**FACS Ca<sup>++</sup> calibration**

1. Prepare solutions of known Ca<sup>++</sup> concentration with Calcium Calibration Buffer Kit #1 (Molecular Probes C-3008):

| Buffer                                 | Add                | to volume of B | Total [Ca <sup>++</sup> ] (mM) | [Ca <sup>++</sup> ] <sub>i</sub> (nM) |
|--|--------------------|----------------|--------------------------------|---------------------------------------|
| 10 mM K <sub>2</sub> EGTA ( <b>A</b> ) |                    |                | 0                              | 0                                     |
| <b>2</b>                               | 4 mL of <b>A</b>   | 1 mL           | 2.00                           | 27                                    |
| <b>3</b>                               | 3.6 mL of <b>2</b> | 1 mL           | 3.74                           | 65                                    |
| <b>4</b>                               | 3.2 mL of <b>3</b> | 1 mL           | 5.23                           | 120                                   |
| <b>5</b>                               | 2.8 mL of <b>4</b> | 1.2 mL         | 6.67                           | 218                                   |
| <b>6</b>                               | 2.4 mL of <b>5</b> | 1.4 mL         | 7.89                           | 410                                   |
| <b>7</b>                               | 2.0 mL of <b>6</b> | 1.8 mL         | 8.89                           | 876                                   |
| 10 mM CaEGTA ( <b>B</b> )              |                    |                | 10                             | ~33 µM                                |

2. Wash Indo-1-loaded cells 1x CLB, and resuspend aliquots of 5 x 10<sup>5</sup> cells in 1 mL each of the buffers.
3. Incubate 90 min 37°C to equilibrate, then collect steady-state indo-1 ratio.

Susannah D. Barbee  
CIT Division of Biology  
barbee@caltech.edu  
3/27/02