

Protocol For CaCl_2 Preparation of Competent E. coli for Transformation with Plasmid DNA

1. Autoclave 1 liter of LB broth in 2 liter baffled flask and cool to room temperature.
2. Inoculate with 5mls-10mls O/N Culture of E. coli started from a single colony.
3. Incubate at 37°C with shaking until the cells are in logarithmic growth 0.4-0.8 at 600nm optical density.
* Should take 4-6 hours.
4. Wash 500mls bottle with 70% EtOH and drip dry inverted and cool on ice.
5. Transfer cells to bottles 2/liter 500mls each
6. Centrifuge at 4000rpm at 4°C keeping cells cold at all times. This is extremely important.
7. Resuspend pellet of cells in 10mls CaCl_2 solⁿ ice cold. Keep on ice for 5 minutes (use 40mls orbe tubes)
8. Centrifuge at 4000rpm 4°C . Resuspend pellet in 10mls CaCl_2 solⁿ ice cold. Keep cells on ice for 30 minutes.
9. Centrifuge at 4000rpm 4°C and drain pellet as dry as possible trying to keep cold.
10. Resuspend pellet in 2.5mls Cold CaCl_2 solⁿ per tube (\approx 500mls cells)
11. Cells should be resuspended well. Keep on ice for a few hours up to 12-24 hours
12. Aliquot in 50 μ l aliquots. 100 aliquots per liter of cells
13. Estimate Competency using plasmid DNA

CaCl_2 Solⁿ :
pH 7.0
8.8g $\text{CaCl}_2 \cdot 2\text{H}_2\text{O}$
3.46g PIPES
150mls Glycerol (100%)
850mls ddH₂O

1000mls

* Filter sterilize 0.2 μ m

* Cells should be aliquoted on ice and frozen immediately after aliquots.