

Protocol for Dideoxy Sequencing Using Sequenase

1. Rxn mixture:
- | | |
|------------------|---------------------|
| DNA | - 4 μ l (26 ng) |
| H ₂ O | - 16 μ l |
| 2M NaOH | - |
| 2M EDTA | - |
| | -> 2 μ l |
| | <hr/> |
| | 22 μ l |

- incubate at 85°C for 5 mins.

2. Put on ice. Add 10 μ l of 0.9M NaOAc
3. Add 1 μ l Primer (200 ng + can dilute 100 μ M (2mer 1/10 use 1 μ l) on the side and mix well
4. Add 95 μ l 100% EtOH. Incubate @ -80°C 15 mins
5. Spin in cold room 20 mins top speed.
6. Rinse with 70% EtOH (100 μ l). Spin 5 mins. Dry in Speed Vac.
7. Resuspend in 10 μ l 1X buffer for Sequenase. mix well.
8. Add:
- | | |
|---------------------|---|
| DTT | - 1 μ l |
| dLab. mix | - 2 μ l (1:5) 1/2 to read close to primer |
| ATP ³⁵ S | - 0.6 μ l |
| Sequenase | - 2 μ l (1:8 in enz. dil. buffer) mix well. |
| | <hr/> |
| | 15 μ l |
- Incubate @ 37°C (R. temp) for 5 mins.

9. Label G, A, T, C tubes (use 0.5 ml Tubes)
10. Add 2.5 μ l of ddNTP to each tube
11. Add 3.5 μ l extension mix to each tube (G, A, T, C) and mix well.
12. Incubate @ 37°C 5 mins. Add 4 μ l stop Solⁿ. Freeze

For SS DNA Sequencing:

| | |
|--------------|-------------|
| template DNA | - 7 μ l |
| 5X Seq. Buff | - 2 μ l |
| -40 primer | - 1 μ l |
| | <hr/> |
| | 10 μ l |

* Heat @ 65°C 2 mins - in 500 μ l beaker then cool down in cold room. Proceed from Step 8.