

# Fluorescence In Situ Hybridization Protocol (for biotinylated probes)

- Denature the chromosomes
  - a. Treat the slide with 70% Formamide, 2XSSC for 2 minutes at 70°C.
  - b. Dehydrate the slide through 70%, 85%, and 100% ethanol for one minute each.
  - c. Air-dry and pre-warm the slide at 37°C on a slide warmer prior to hybridization.
- Denature the Probe
  - a. Mix 7 ul of hybridization buffer with 3 ul probe mixture in a 0.5 ml tube.
  - b. Denature the hybridization mix at 70°C for 5-10 minutes.
- Hybridization
  - a. Put the hybridization mix onto the slide and apply a 22x22mm coverslip.
  - b. Seal the edge of the coverslip with rubber cement.
  - c. Keep the slide in a moist chamber overnight at 37°C.
- Fluorescence staining
  - a. Peel off the rubber cement and take off the coverslip.
  - b. Wash the slide with 50% formamide, 2X SSC at 45°C three times for 5 minutes.
  - c. Wash the slide with 4X SSC, 0.05% Tween20 three times for 2 minutes.
  - d. Wash the slide for 2 minutes in 4X SSC at room temperature.
  - e. Put 100 ul of Avidin (5 ug/ml Avidin in PNM buffer) onto the slide and apply a coverslip.

- f. Incubate the slide for 20 minutes at room temperature in a DARK moisture chamber.
  - g. Repeat steps (c) and (d)
  - h. Put 100 ul anti-Avidin (1:200 dilution in PNM buffer) and apply coverslip.
  - i. Incubate the slide for 20 minutes at room temperature in a DARK moisture chamber.
  - j. Repeat steps (c), (d), and (e), followed by step (c) and (d).
  - k. Put 40 ul of antifade solution (0.5 ug/ml PI or DAPI as counterstain) onto the slide and apply coverslip.
- Examine the slides or store in the Dark

Home