

Gel Shift (EMSA) Conditions

Gel Shift Recipe

| | | |
|---------------------|-------------|-----------------------------|
| H ₂ O | - 34.5ul | - incubate at R.T. 15 mins. |
| G.S | - 2.5ul | - add 20ul STE buffer |
| (300ng) DNA | - 3.0ul | - NucTrap |
| NTP's | - 5.0ul | - add 12ul STE buffer |
| ³² PdATP | - 4.0ul | Final Vol. = 240ul |
| Klenow | - 1.0ul | |
| | <u>50ul</u> | - Current 1ul probe |

EMSA Buffer Conditions

Crude Nuclear Extract Buffers:

Protease Inhibitors:

| | | |
|-------------|----------|---|
| pepstatin A | 0.5ug/ml | 50 mM KCl |
| leupeptin | 5ug/ml | 4 mM MgCl ₂ |
| Chymostatin | 5ug/ml | 20 mM K ₂ PO ₄ pH 7.4 |
| antipain | 5ug/ml | 1 mM β -mercaptoethanol |
| aprotinin | 5ug/ml | 20% glycerol |
| benzamidine | 5 mM | |

- dialyse extract at 4°C for 30 mins in 500x volumes of extract buffer shown above.
- Store in small aliquot at -70°C. Use 5ug per rxn if possible.

Make 4x EMSA buffer stock:

| | | | |
|--------------------|---|----------------------|-----------------------------------|
| 40 mM HEPES pH 7.4 | } | To 500ul add | - 10ul BSA (20mg/ml) |
| 20% glycerol | | - 2.5ul PMSF (100mM) | |
| 4 mM EDTA | | - 2.0ul DTT (1.0mM) | |
| | | To each rxn | { - 2.0ul dIdc (1ug/ml) = 2ug/rxn |
| | | 50ul | { - 2.0ul KCl (1M) = 40mM |

- Add components of each rxn to tube 1.5ul eppendorf and incubate at room temperature for 5 minutes
- add probe (20ng) and incubate further at R. temp for 15 minutes.
- Run gel 160 volts ... Perm 165V 16w or 50V o/p

* Varying the concentration of non-specific competitor dIdc from 20ng to 2ug per rxn can allow the interaction of specific protein-DNA complexes.