

# Gelatin Zymography Protocol

## Key Steps:

- 1) Transfect MMP producing cell line (C33a) with expression vector containing protein of interest.
- 2) Collect media conditioned with MMPs
- 3) If desired concentrate protein in media
- 4) Run on 10% SDS-PAGE with 0.1% gelatin
- 5) Incubate gel for 24 hours and then stain

## Day 1

- 1) Plate C33a cells in appropriate size plate so that they will be 80% confluent the following day. (Works best using 10cm plates)

## Day 2

- 1) Transfect C33a cells when 80 percent confluent using Lipofectamine 2000 and incubate for 24 hours.
  - \* Only use media w/o FBS and Pen/Strep
    - Using Pen/Strep with Lipo 2000 will kill the cells
    - FBS contains Bovine MMPs which will show up on the zymograph

## Day 3

- 1) Make a 10% SDS-PAGE gel containing 0.1% (1mg/ml) gelatin in the resolving gel.
  - \* Do not use buffers containing EDTA.
  - \* Make stock gelatin fresh every time at a concentration of 10mg/ml.
  - \* Use gelatin immediately after making.
- 2) Collect the media off of the cells and spin at 10,000 rpm for 5mins to remove any cell debris.
- 3) Concentrate approx. 6mls of media down to 300ul using Centricon 10.
- 4) Mix 20ul of media with 2x SDS loading buffer w/o 2-mercaptoethanol.
  - \* Optional: incubate the samples at 37°C for 20 mins.  
Might help activate the MMPs???
- 5) Run gel at 4°C at 150-200V until the 50kD marker leaves the gel.
- 6) Cut off marker lanes
- 7) Wash gels 4 times for 30 mins each at room temperature in washing buffer.
  - \* The buffer removes the SDS from the gel allowing the protein to refold
- 8) Wash gels 30mins in ddH<sub>2</sub>O
- 9) Incubate the gel for at least o/n (generally 24-48h) at 37°C in the incubation buffer.

## Day 4

- 1) Quickly rinse the gel with distilled water.
- 2) Stain the gel with Comassie Brilliant Blue solution for 1 hour.
- 3) Destain with destain solution for 30mins-1hr or until there is good resolution between the bands and blue background.

4) Image the gel using the Fuji Imager BEFORE drying.

## **Buffers and Reagents**

### **Stock Gelatin**

0.1 g gelatin  
ddH<sub>2</sub>O to 10ml

### **5X Resolving Gel Buffer pH 8.8**

63.58g Tris  
2.5g SDS  
ddH<sub>2</sub>O to 250ml  
Adjust pH to 8.8

### **5X Stacking Gel Buffer pH 6.8**

47.75g Tris  
6.25g SDS  
ddH<sub>2</sub>O to 250ml  
Adjust pH to 6.8

### **SDS-Loading Buffer**

1.25ml 1M Tris, pH 6.8  
2.0ml 20% SDS  
3.75ml ddH<sub>2</sub>O  
2.0ml Glycerol

### **Washing Buffer**

25ml 1M Tris, pH 7.6  
5ml 1M CaCl<sub>2</sub>  
0.2g NaN<sub>3</sub>  
500ul 0.5M ZnCl<sub>2</sub>  
12.5ml Triton X  
ddH<sub>2</sub>O to 1L

### **Incubation Buffer**

25ml 1M Tris, pH 7.6  
5ml 1M CaCl<sub>2</sub>  
0.2g NaN<sub>3</sub>  
75ml 5M NaCl  
ddH<sub>2</sub>O to 1L

### **Comassie Brilliant Blue**

40% MeOH  
10% Acetic Acid  
0.1% CCB

### **Destain Solution**

40% MeOH  
10% Acetic Acid