

Page Six of Six
Catalog # 17-281
Lot # 18895

HDAC Protocol (³H labeled Histone H4 Peptide)

Stock Solutions:

1. **³H labeled Histone H4 Peptide:** Prepared according to above protocol; 20,000 CPM per assay is recommended.
2. **5X HDAC Assay buffer:** Add PMSF to a final concentration of 5mM prior to assay.
3. **Quenching Solution:** Add 259 μ l of HCl and 28 μ l of acetic acid to 2713 μ l of distilled water. Prepare immediately before use.

Assay Protocol

1. Label microcentrifuge tubes (in duplicate or as appropriate) for samples and controls (see Technical Note a).
 2. To each tube add :
 - 40 μ l of 5X HDAC Assay Buffer stock solution with PMSF
 - 20,000 CPM [³H]-acetyl Histone H 4 peptide
 - Source of HDAC (see Technical Note b)
 - distilled water to 200 μ l
- Prepare two control reactions. Add 10 μ g HeLa Nuclear Extract to each tube. In one tube add 10-50 μ l 1M Sodium Butyrate (see Technical Note c)
3. Centrifuge briefly to collect the components in the bottom of the tube and incubate reactions on a rocker at room temperature for up to 24 hours, or at 37°C for several hours.
 4. Centrifuge briefly to collect the components in the bottom of the tube.
 5. Add 50 μ l Quenching Solution to stop the reaction and vortex.
 6. In a chemical fume hood, extract released [³H]-acetate. Add 600 μ l of ethyl acetate to each tube. Vortex and centrifuge samples for one minute at 14,000 x g in a microfuge to separate phases.
 7. Determine radioactivity in two 200 μ l aliquots of each ethyl acetate phase. Transfer each aliquot to a separate scintillation vial containing scintillation fluid, mix thoroughly and measure CPM. Compare CPM between samples incubated with or without sodium butyrate.

Technical Notes:

- a. For each test sample, a second duplicate assay should be performed in the presence of 50-250mM Sodium Butyrate to demonstrate specificity of deacetylation.
- b. HDAC activity is present in assayable levels in nuclear extracts and in cell lysates prepared with RIPA buffer. HDAC-containing immune complexes, washed in 1X HDAC Assay buffer may also be assayed with this protocol. Use 5 μ g anti-HDAC1 (Catalog # 06-720) or 10 μ g anti-HDAC3 (Catalog # 06-890) to immunoprecipitate the respective HDACs from 500 μ g of lysate or nuclear extract.
- c. For assays performed with HeLa nuclear extract, 50mM Sodium Butyrate inhibited 70-75 % of the HDAC activity.