

## Hoechst 33342 Staining of Viable Cells for cell cycle analysis/sorting

---

- For cell cycle analysis and sorting of unfixed samples
  - Requires ultraviolet argon laser excitation, 450nm emission
  - Binds preferentially to A-T base regions in DNA
  - Non-intercalative DNA binding
  - DiOC5 can increase resolution of DNA distribution
- 

### REAGENT LIST:

- Hoechst 33342 Staining solution (HO)  
Hoechst 33342 1.0 mg/ml in dH2O
  - DiOC5 (Molecular Probes) 1.0mg/ml in DMSO
- 

### PROTOCOL:

#### Staining:

- 1.HO is added to cells in culture medium at from 1.0-5.0 ug/ml.
  - 2.Cells are incubated in HO at 37 degrees C for 30-60 minutes.
  - 3.Cells are analysed without washing while in the media containing the HO.
- 

### TIPS:

- 1.DO NOT SUBSTITUTE HOECHST 33258 FOR VIABLE CELL STAINING.
  - 2.HO stock should not be mixed in PBS, the dye will precipitate.
  - 3.Stock dye solution may be frozen and is stable for at least one month in the refrigerator in foil.
  - 4.Adherent cells should be stained in vitro, then trypsin and trypsin-neutralizing solutions used should contain the same dye concentration
  - 5.Improved resolution (CV'S) can be obtained by adding 0.1-0.3ug/ml DiOC5 at the same time as HO incubation.
- 

[Back to the DFCI Core Flow Cytometry Home Page](#) •