

*Rajeev,  
Here is the processivity protocol.  
Please place copy in protocol book  
for lab. T.S. Ede*

Ruby Chan  
06-12-06

**In vitro Processivity Assay**

1. end-label d(T)<sub>16</sub> with [ $\gamma$ -<sup>32</sup>P]-ATP (X 2 tubes)
  - 30 $\mu$ l oligo d(T)<sub>16</sub> (0.1 $\mu$ g/ $\mu$ l stock, 3 $\mu$ g)
  - 5 $\mu$ l 10X T4 PNK buffer
  - 10 $\mu$ l [ $\gamma$ -<sup>32</sup>P]-ATP (3000Ci/mmol, 10mCi/ml)
  - 1 $\mu$ l 10U/ $\mu$ l T4 PNK
  - 4 $\mu$ l ddH<sub>2</sub>O
  - 50 $\mu$ l total
2. incubate at 37<sup>o</sup>C for 30min
3. stop reaction by adding 2 $\mu$ l 0.5M EDTA
4. transfer 52 $\mu$ l into a Bio-Rad P-6 micro spin-column, cf 4min
5. pool 2 flow-through together, store at -20<sup>o</sup>C
6. mix 25 $\mu$ l poly(dA) (25 $\mu$ g) with 50 $\mu$ l labeled oligod(T)<sub>16</sub>
7. heat at 70<sup>o</sup>C, 5min, xfer to RT, cool for 1/2hr
8. pass annealed poly(dA) template thru 2 P-30 micro spin-columns (aliquot template into 2 columns so that centrifuge is balanced)
9. set up TNT samples and incubate on ice for 30min

*732-6221  
micro Bio-spin 6 columns*

sample

Pol Assay Soln.

	<u>stock</u>	(10X) <sup>-</sup> <u>41 <math>\mu</math>l</u>	
100mM (NH <sub>4</sub> ) <sub>2</sub> SO <sub>4</sub>	2M	12.5 $\mu$ l	
20mM Tris-HCl (pH7.5)	1M	5	
3mM MgCl <sub>2</sub>	100mM	7.5	
0.1mM EDTA	10mM	2.5	
0.5mM DTT	50mM	2.5	
4% glycerol	100% <i>sc</i>	10	20
40 $\mu$ g/ $\mu$ l BSA	10 $\mu$ g/ $\mu$ l	1 $\mu$ l	
		<del>50</del>	
		<u>51 <math>\mu</math>l</u>	(5X)

10. add 3 $\mu$ l poly(dA):oligo(dT) and incubate for 10min at RT
11. add 15 $\mu$ l Pol Assay Soln. (4.1 $\mu$ l) + 1mM dTTP(2.5 $\mu$ l) (final conc.=100 $\mu$ M) + ddH<sub>2</sub>O (8.4 $\mu$ l)
12. incubate for 1hr, 37<sup>o</sup>C

13. add 25 $\mu$ l ddH<sub>2</sub>O and 50 $\mu$ l phenol:chloroform:isoamylalcohol (25:24:1), mix, cf, xfer aqueous phase to new tube
14. extract with 50 $\mu$ l chloroform:isoamylalcohol (24:1)
15. add 20 $\mu$ l urea-PAGE loading dye
16. pre-run 15% Urea-PAGE at 300V for 30min
17. transfer 10 $\mu$ l sample into new tubes
18. flush wells with 1X TBE to get rid of urea
19. while flushing the wells, heat sample (10 $\mu$ l) at 98<sup>o</sup>C for 5min, immediately load onto gel
20. resolve samples in 15% Urea-PAGE at 300V for 2hrs and 45min (regular size) or at 500V for 8hrs (long gel)

15% Urea-PAGE (regular size)

	<u>40ml</u>
40% acry/bis (19:1)	15ml
7M urea	16.8g
10X TBE	4ml
ddH <sub>2</sub> O	7.77ml
10% APS	300 $\mu$ l
TEMED	30 $\mu$ l

- dissolve urea in acry/bis, TBE, and ddH<sub>2</sub>O at 42<sup>o</sup>C
- cool gel soln at RT for 0.5 hr, add TEMED and APS last
- gel starts polymerizing after 5min, and let gel set for at least 15min

15% Urea-PAGE (61 cm long)

	<u>80ml</u>
40% acry/bis (19:1)	30ml
7M urea	33.6g
10X TBE	8ml
ddH <sub>2</sub> O	15.65ml
10% APS	500 $\mu$ l
TEMED	50 $\mu$ l

- dissolve urea in acry/bis, TBE, and ddH<sub>2</sub>O at 42<sup>o</sup>C
- cool gel soln at RT for 0.5 hr, add TEMED and APS last
- pour gel soln. with 25ml pipet in 5min, gel starts polymerizing after 5min
- let gel set for at least 15min
  
- if not dry gel, put gel on filter paper and wrap gel with saran wrap
- bands start diffusing after 1 day if gel is not dried
- for long gel, expose overnight with blue film (Fuji film)
- for regular-size gel, expose for 2hrs with Kodak X-omat

**Ruby Chan**  
**06-12-06**

denaturing PAGE loading dye (10ml)

9ml deionized formamide  
2ml 0.5M EDTA (pH8) (final=10mM)  
0.1g bromophenol blue (0.1% final)  
0.1g xylene cyanol (0.1% final)

10X TBE

890mM Tris base	108g
890mM boric acid	55g
20mM EDTA (stock=0.5M)	40ml
add ddH <sub>2</sub> O to 1L	




Fig.: Processivity assay showing the efficiency of replication machinery of BC3 and BJAB cells. The nuclear extracts from these cells were incubated with poly(dA):oligo(dT)<sub>16</sub> for 30 minutes (lanes 2 and 3 respectively). DNA products were fractionated by 7M urea-15% PAGE. Poly(dA):oligo(dT)<sub>16</sub> template was electrophoresed in lane 1 to show the size of primers (16 base). Nonspecific DNA products are labelled by closed asterik symbol and specific are marked by open asterik symbol.