

### ***IP/HAT assay***

Lyse cells in  $\sim$ 0.5% NP-40 lysis buffer

Preclear for 2 hours with protein A agarose

Incubate with antibody overnight at 4°C

Add protein A agarose and incubate for an additional 2 hours at 4°C

Wash immunoprecipitates 2 times with lysis buffer and 2 times with HAT assay buffer

Incubate beads in 30  $\mu$ l HAT assay buffer, histones and [ $^3$ H]acetyl CoA at 30°C for 1 hr

*HAT assay buffer:*

50 mM Tris-HCl (pH 8.0) ✓

10 mM sodium butyrate ✓

10% glycerol ✓

1 mM DTT ✓

1 mM PMSF ✓

+ 25  $\mu$ g crude histones *no type II A*

+ 50 nCi [ $^3$ H]acetyl-CoA

*3.2 Ci/mmol*

Spot 20  $\mu$ l of reaction onto Whatman P-81 filter paper

Wash for 30 min with 0.2 M sodium carbonate buffer (pH 9.2) (Change buffer 2 times)

Dry filters and count in liquid scintillation counter