

Isolation of mRNA from Suspension of Cells to RNazol B

1. Pellet cells: $100 \times 10^6 - 1 \times 10^7$ cells into a tight pellet at 4°C at 1500rpm for 8-10 mins.
* cells can be frozen here at $-80^\circ\text{C} - -140^\circ\text{C}$ before proceeding.
 2. Resuspend cells thoroughly on ice in 1ml ice cold PBS
 3. Add RNazol B 20ul / 100 mill - 500 mill (pipette up and down few times) and cover tightly. Shake vigorously for 15-30 Secs (do NOT vortex) and let them stay on ice or 4°C for 5 mins after addition of 2ul / 20ul lysate of CHCl_3
 4. Centrifuge at 4°C 12,000 xg for 15 mins (9000rpm in SA-600: use 7500rpm)
 - * 5. For larger amount of cells remove bottom organic phase and repeat step 3 by adding 20ul of RNazol B and passing thru pipette a few times. Add 2ul CHCl_3 then shake vigorously for 15 Secs. let sit on ice/ 4°C 5 mins.
 6. Centrifuge at 4°C 12,000 xg for 15 mins (7,500rpm in SA-600)
 7. Remove aqueous phase and transfer to a fresh tube on ice and add an equal volume of Isoopropanol and store samples at 4°C for 15 mins / or overnight here.
 8. Centrifuge samples for 15 mins at 7,500rpm at 4°C .
 9. Remove supernatant and wash pellet once in 75% EtOH by vortexing and centrifugation at 7500rpm 4°C . Use 10ul 75% EtOH.
 10. Dry pellet briefly (* DO NOT OVERDRY) by inversion on bench top.
 11. Dissolve the RNA pellet in 0.5% SDS in 1ml EDTA pH7. by vortexing and passing thru pipette a few times.
 12. Incubate at 65°C for 10-15 mins
 13. ppt. RNA in 0.2M NaCl to one vol. isopropanol / or 2 vols. EtOH for 1 hour at -20°C .
 14. Wash pellet again in ^{10ul} 75% EtOH: Vortex & Centrifuge at 4°C , 7,500rpm for 15 mins.
 15. Resuspend pellet in appropriate vol. of RNase free TE (H₂O DEPC treated)
 16. Take Spec. Reading for Concentration.
* For polyA prep.
- * Aliquot into separate tubes if not using immediately.
Store at -80°C