

# Isolation of Total RNA from Cultured Cells

## I. Solutions

*make up fresh each time*

### A. GIT buffer

- 9.452 g  
.334 ml  
19.66 ml  
.334 ml
1. 47.26g guanidine isothiocyanate (4M)
  2. 1.67ml 3M sodium acetate, pH 6.0 (50mM)
  3. dH<sub>2</sub>O to 98.33ml
  4. filter through 0.2um filter
  5. 1.67ml 2-mercaptoethanol
  6. store in the dark at RT

### B. CsCl buffer

1. 47.98g CsCl (5.7M)
2. 0.415ml 3M sodium acetate, pH 6.0 (25mM)
3. dH<sub>2</sub>O to 50ml
4. filter through 0.2um filter
5. store at RT

### C. 2X PK buffer

1. 1.21 g Tris (0.2M); 10ml 1M Tris, pH 7.5
2. 0.46g EDTA (25mM); 2.5ml 0.5M EDTA
3. 0.88g NaCl (0.3M); 3.75ml 4M NaCl
4. 1.0g SDS (2%)
5. dH<sub>2</sub>O to 45ml
6. pH to 7.5 with HCl
7. ~~1ml 10mg/ml proteinase K (200ug/ml)~~
8. dH<sub>2</sub>O to 50ml
9. ~~store at -20°~~

*Use proteinase K  
at 100 µg/ml*

## II. Procedure

- A. Pellet 400-500ml of cells (10<sup>5</sup>-10<sup>6</sup> cells/ml) in the Sorvall T600B centrifuge, 5 min., 4°, 1000rpm.
- B. Decant supernatant and wash cell pellets with 10ml cold (4°) Dulbecco's PBS. Pellet cells as above.
- C. Resuspend cells in 18ml GIT buffer and vortex vigorously.
- D. Shear the chromosomal DNA by transferring the cells to a 60ml syringe and shearing through a 20 gauge needle (3X).
- E. Into 2 SW41 centrifuge tubes add 3ml CsCl buffer and layer 9ml of the GIT-cell mixture. Pellet the RNA by centrifuging 21 hr., 32K rpm, 20°. *The temperature is critical since the CsCl will precipitate at lower temperatures.*
- F. Remove the GIT buffer layer and part of the CsCl layer using a Pasteur pipette and house vacuum. Decant the remaining CsCl layer using a Pasteur pipette being careful not to disturb the RNA pellet. The pellet should appear clear.
- G. ~~Rinse the pellet with 1ml cold (-20°) 80% EtOH by swirling gently.~~ Do not rinse the portion of the tube that contained the GIT buffer. Carefully remove the EtOH with a Pasteur pipette.

*use 2N NaOH to wash & flush well w/ dH<sub>2</sub>O*

*Andrew has used brake (no setting) maybe better to use any brake (zero setting)*

*Dry the pellet + resuspend in TE + .1% SDS  
Add NaOAc + EDTA  
(rinse, dry) then resuspend in 2X PK bu  
originally from pellet*

*ddH<sub>2</sub>O + .1% SDS to resuspend originally from pellet*

*TE.*

- H. Air dry the pellet and then add 1ml of 1X PK buffer and incubate at 50° for 30-60 min. to completely solubilize the RNA.
- I. Add an equal volume of buffered phenol (see Maniatis p 438), mix by vortexing and separate the phases by spinning in the microfuge, 2 min, RT, full speed.
- J. Transfer the upper aqueous phase with a P-200 and re-extract 2-3X until the interphase is clear.
- K. Add 1/10 volume 3M sodium acetate, pH 6.0, and 2.5 volumes cold (-20°) EtOH. Mix by inverting several times and place on dry ice for 20 min. Pellet the RNA by spinning in a microcentrifuge, 15 min, 4°, full speed.
- L. Decant the EtOH and dry the pellet in the speed-vac for 5 min., ~~heat on~~. Resuspend the dried pellet in 200-500ul dH<sub>2</sub>O and store at -70°. For long term storage, pellet the RNA as above and store as is at -70°.

transfer to eppendorf

then split into 2 eppendorf  
to do phenol

### III. References

- A. Basic Methods in Molecular Biology, Section 11-1, Davis, L.G., Dibner, M.D., and Battey, J.F. (1986).
- B. Molecular Cloning, A Laboratory Manual, Chapter 6, p. 196, Maniatis, T., Fritsch, E.F., and Sambrook, J. (1982)
- C. Current Protocols in Molecular Biology, Vol. 1, Section 4.2.1-4.2.5

Tissues

- see notes
- (resuspend in 100) d<sub>2</sub>O