

Lambda Packaging in PLK-A

- from a single colony inoculate 5mls of LB supplemented w/ $MgSO_4$ and Maltose and grow o/n at $37^\circ C$ w/ shaking.
- dilute culture 1/10 in LB w/ Maltose in 5mls and grow for 4-6 hours or until density at 600nm $\approx 0.4-0.6$ O.D.
- pellet bacteria at 2000 rpm 10 mins. Table top 10mins R.T.
- Resuspend cells in half original volume 10mM $MgSO_4$.
- * Cells here can be stored 2-3 days before use

Packaging Protocol:

- 1- Remove appropriate number of packaging tubes (Red & Yellow) or (blue and white) from Stratagene GigaPak IXL extracts and place on ice immediately.
- 2- Aliquot DNA 1-4ul (0.1-5ug DNA) in tube. Add contents of 10ul (red/blue) tube to DNA tube and pipette up and down to mix and place on ice.
- 3- Quickly add 15ul of (yellow/white) tube to tube w/ mixture and pipette up and down to mix. Microfuge briefly to get rid of bubbles (5 secs)
- 4- Incubate tube at R.T. $22^\circ C$ for 2 hours
* Do NOT EXCEED 2 HOURS (9 mins optimum).
- 5- Add 500ul SM Buffer and 10ul $CHCl_3$. Mix contents of tube gently. Centrifuge to sediment debris briefly (30secs).
- 6- Store supernatant w/ virus at $4^\circ C$ for use.

Infection:

- Bacteria (PLK-A) resuspend in $MgSO_4$ 10mM 200ul - 25ul phage
200ul - 50ul phage (up to 100ul of phage)
- mix contents in each tube and incubate at $37^\circ C$ for 20 mins.
- Add 1ml LB w/ maltose. Shake at $37^\circ C$ for 1 hour.
- plate on ampicillin plates. Incubate o/n $37^\circ C$

Lambda Packaging in PLK-A

- From a single colony inoculate 5mls of LB supplemented in MgSO₄ and Maltose and grow overnight at 37°C w/ shaking.
- Dilute culture 1/10 in LB w/ Mg + Maltose into 5mls and grow for 4-6 hours or until density at 600nm = 0.4-0.6 O.D.
- Pellet bacteria at 2000rpm 10 mins. Table top for 10 mins R.T.
- Resuspend cells in half original volume 10mM MgSO₄
- * Can be stored for 2-3 days before use.

Packaging Protocol

1. Remove appropriate number of packaging tubes Red & yellow or blue & white for Gigapack II XL extracts and place on ice immediately
2. Aliquot DNA 1-4 μl (0.1-5 μg DNA) in tube. Add contents of (10 μl) Red (or blue) tube to DNA and pipette up & down to mix and place on ice
3. Quickly add 15 μl of yellow (white) tube to DNA tube and pipette up and down to mix. Vortex briefly to get rid of bubbles (3-5 secs)
4. Incubate tube at R.T. 22°C for 2 hours * Do NOT EXCEED 2 HOURS. 90 mins optimum.
5. Add 500 μl SM Buffer and 10 μl clets. Mix contents of tube gently. Centrifuge to sediment debris.
6. Store Supernatant at 4°C for use.

Infection

1. Bacteria Resuspend in MgSO₄ 10mM

	<u>Bacteria</u>	<u>Phage</u>
-	200 μl	25 μl
-	200 μl	50 μl
-	200 μl	75 μl
-	200 μl	100 μl

- Mix contents in each tube.
2. Incubate at 37°C for 20 minutes. Add 1 ml of LB w/ Mg + Maltose medium. Shake at 37°C for 1 hour
3. Plate on LB ampicillin plates. Incubate O/N 37°C

① Recipe For Antifade

pH to 8.0

1x PBS - 4 mls

Phenylenediamine HCl - 40mg

Glycerol - 36 mls

1M bicarb/carb - 1.5 mls pH 9.5

② MgSO₄ for Phage Packaging

10 mls - 1M MgSO₄

90 mls - ddH₂O

Filter sterilize. Volume 100 mls.

③ LB w MgSO₄ and Maltose

10g - Bacto tryptone

5g - yeast Extract

5g - NaCl

1000 mls - ddH₂O autoclave for 20 minutes at 121°C 15 lbs.

10 mls - 1M MgSO₄

20 mls - 10% maltose

} add after cooling by filter sterilizing thru 0.2um Filter.

④ SM Buffer

5.8g - NaCl

2.0g - MgSO₄·7H₂O

50.0 mls - 1M Tris. pH 7.5

5.0 mls - 2% w/v Gelatin

↑ 1 liter - dd H₂O autoclave