

Nuclear Protein Extract

6/30/95

Buffer A

10mM HEPES pH 7.9 - 100ul 1M
 10mM KCl - 100ul 1M
 1.5mM MgCl₂ - 15ul 1M
 dH₂O - to 10mls
10mls

Buffer B

20mM HEPES pH 7.9 - 200ul 1M
 10% glycerol - 1000ul 100%
 420mM NaCl - 840ul 5M
 1.5mM MgCl₂ - 15ul 1M
 0.2mM EDTA - 4ul 0.5M
 dH₂O - to 10mls
10mls

Add 1M DTT - 50ul/10ml
 100mM PMSF - 50ul/10ml
 10mg/ml Aprotinin - 10ul/10ml

- add jnt before use to both A and B buffers on ice
- Collect cells ~ 50 → 100 x 10⁶ at 1-2000 rpm 10 mins 4°C
- Wash cells once in 10-25 mls PBS. Pellet as above
- resuspend in 4 vols. to size of pellet Buffer A ≈ 500-1000ul.
- Incubate on ice for 1 hour.
- Transfer to douncer and homogenize 20-30 strokes
- transfer to 1.5 ml epp. tube. Centrifuge @ 2000 rpm 5 mins 4°C
- Aspirate and resuspend in 1 ml Buffer A.
- centrifuge @ 2000 rpm 5 mins, 4°C. Aspirate.
- Resuspend in 3 vols. pellet with buffer B.
- incubate on ice 30 mins.
- Microfuge 20 mins @ 13,000 rpm 4°C.
- Transfer to fresh 1.5 ml epp.
- Snap freeze. -80°C

Bradford Assay:

Dilution 1:200

{	795ul H ₂ O
	5ul protein (NE)
	200ul Bradford
	<u>1000ul</u>

Buffer A Extract Buffer

100 mM	Tris-cl	pH 9.0	<u>10 ml</u> 1 ml	1 mM	Tris-cl	pH 9.0
100 mM	NaCl		1 ml	1 mM	NaCl	
5 mM	KCl		50 μ l	1 mM	KCl	
0.5 mM	MgCl ₂		5 μ l	1 mM	MgCl ₂	
1 mM	CaCl ₂		10 μ l	1 mM	CaCl ₂	
0.5%	NP40		500 μ l	10%	NP40	

Increase to 10 ml with ddH₂O

Buffer B Gradient Buffer

10 mM	HEPES	pH 7.8	<u>50 ml</u> 500 μ l	1 mM	HEPES	pH 7.8
5 mM	Na ₂ HPO ₄		250 μ l	1 mM	Na ₂ HPO ₄	
5 mM	KCl		250 μ l	1 mM	KCl	
0.5 mM	MgCl ₂		25 μ l	1 mM	MgCl ₂	
1 mM	DTT		500 μ l	100 mM	DTT	
0.5%	NP40		(250 μ l	0.5%	NP40	(final conc)