

Nuclear Protein Extraction

7/7/96

Buffer A

10mM HEPES pH 7.9 - 100ul 1M
 10mM KCL - 100ul 1M
 1.5mM MgCl₂ - 15ul 1M
 dH₂O sterile - to 10mls

Buffer B

20mM HEPES pH 7.9 - 20ul 1M
 10% glycerol - 1ml 100%
 420mM NaCl - 840ul 5M
 1.5mM MgCl₂ - 15ul 1M
 0.2mM EDTA - 4ul 0.5M
 dH₂O sterile - to 10mls

Buffer C

20mM HEPES pH 7.9 - 20ul 1M
 30% glycerol - 3ml 100%
 1.5mM MgCl₂ - 15ul 1M
 0.2mM EDTA - 4ul 0.5M
 dH₂O sterile - to 10mls

Add: 1M DTT - 50ul / 10ml
 100mM PMSF - 50ul / 10ml
 10mg/ml Aprotinin - 10ul / 10ml

} - to Buffer A, B and C just before use.
 } - place buffers on ice for extended period of use.

- collect 100 million cells at 1-2000rpm 10 mins at 4°C.
- Wash cells once in cold PBS 25mls. Pellet as above.
- Resuspend in 4 vols. to size of the pellet in buffer A = 500ul - 1000ul
- Incubate on ice for 1 hour. Transfer to douncer and homogenize 20-30 strokes
- Transfer to 1.5mls eppendorf tube. Centrifuge @ 2000rpm 5 mins at 4°C.
- Aspirate and resuspend in 1ul Buffer A. Centrifuge @ 2000rpm 5 mins at 4°C and aspirate
- Resuspend in 2 volumes pellet with buffer B.
- incubate on ice for 30 mins. Microfuge 20 mins @ 13,000rpm 4°C.
- transfer to fresh 1.5ml eppendorf. Add equal volume of Buffer C.
- Snap freeze in dry ice ethanol bath. Store at -80°C

Bradford Assay:

795ul dH ₂ O	} Dilution 1:200
5ul protein (NE)	
200ul Bradford	
<u>1000ul</u>	

* Need concentration of Nuclear Extract 3.5mg / ul

Binding Buffer

20 mM phosphate buffer pH 6
10 mM $MgCl_2$
0.1 mM EDTA
2 mM DTT
0.01% NP40
0.1 mM PMSF
50-100 mM NaCl
15% glycerol
100 μ g/ml BSA

* Binding rxns were done in
rxn of binding buffer

- For gel retardation protein fractions
were added to binding buffer
containing DNA labeled probe
(1000-3000 cpm) and poly dI-dC
as non-specific competitor

Protein Samples: Buffer C

20 mM potassium phosphate buffer pH 6.0
10 mM $MgCl_2$
0.1 mM EDTA
1 mM DTT
0.01% NP40
0.1 mM PMSF
10% glycerol
100 mM NaCl

} protein samples dialyzed and
stored in this buffer

- 0.25 μ g 5'-end labeled oligo in 1 μ l rxn vol. w 200 μ g poly dI-dC
and 3-6 μ g of nuclear extract

- with purified protein no poly(dI-dC) but 100 μ g/ml BSA.

* Incubations :
1 - 15 min R.T
2 - 30 min R.T
3 - 1 hour 4°C w antibody before EMSA

For Gel Shift Assay:

Protocol:

Binding Buffer:

20 mM Phosphate buffer pH 6.0,
10 mM $MgCl_2$
0.1 mM EDTA Na_2
2 mM DTT
0.01% PMSF
100 mM NaCl
15% Glycerol
100 μ g/ml BSA

Dialysis Buffer for purified proteins Samples: Buffer C.

20 mM K PO₄ buffer pH 6.0
10 mM $MgCl_2$
0.1 mM EDTA
1 mM DTT * proteins also
0.01% NP40 stored in
0.1 mM PMSF this buffer
10% glycerol
100 mM NaCl

For Gel Retardation Assay:

1. Do rxns in total volume of 20 μ l of binding buffer
2. Amount of probe equivalent to 0.25 μ g 5'-end labeled oligo.
with 200 μ g poly dI-dC in 10 μ l buffer.
3. Amount of protein in total rxn 3-6 μ g
4. When using purified protein no poly dI-dC but 100 μ g/ml BSA
5. Incubate rxns for 30 mins at R-Temp
6. For rxns which includes antibodies incubate for 1 hour at 4°C prior to running EMSA gel.