

# ONE COLOR FISH

## Nick Translation (BioNick Labeling System GibcoBRL)

\*IMPT: Always keep the enzyme in the freezer  
Keep the total reaction volume at 50ul  
[BAC] = 1ug DNA

- |                          |                        |                      |
|--------------------------|------------------------|----------------------|
| 1. Labeling Rxn Example: | BAC DNA                | 5.5ul                |
|                          | 10x Enzyme             | 5.5ul                |
|                          | 10x <del>M</del> dNTPs | 5.0ul                |
|                          | water                  | to 50ul total volume |

Pulse spin labeling reaction  
Nick translate 1 hour at 16°C

2. Check 10ul reaction on a 2.0% agarose gel
  - should be a smear between 200 and 500bp
  - keep nick translation reaction on ice until adding 5ul stop Buffer to stop the reaction
3. Add 5ul Stop Buffer
4. Add 20ug COT-1 DNA and vortex to mix
5. Add 1/10 volume 3M NaOAc and 2 volums 100% EtOH to reaction and mix. Place in a -80°C freezer for 30min.
6. Centrifuge in a 4°C microcentrifuge for 30min. at high speed.
7. Decant supernatant and wash pellet with ice cold 70% EtOH.
8. Dry pellet in Speed Vac for 5-10min.
9. Resuspend pellet in 3ul ddH<sub>2</sub>O and 7ul hybe mix containing:

50% Formamide	3.5	X3
2X SSC	.7	10
10% Dextran Sulfate	2.77	2.1

  - a. Pulse spin
  - b. Vortex to mix and resuspend
  - c. Pulse spin
  - d. Place in 37°C water bath to resuspend pellet (about 10min.)
10. Denature probe by heating to 70°C for 7.5-8min.
11. Preanneal probe at 37°C for the appropriate amount of time:

Cosmids:	45-60min
BAC:	1-5h (2.5-3h)
YAC:	1-5h (2.5-3h)