

Protocol : Purification of *Taq* DNA polymerase

Reference: Pluthero, F.G. *Nucleic Acids Research* (1993) **21**, 4850-4851

Buffers:		per 1 L	
Buffer A	50 mM Tris-HCl pH 7.9	50 ml 1 M pH 7.9	150 ml
	50 mM dextrose	50 ml 1 M	5 ml
	1 mM EDTA	2 ml 0.5 M	2 ml
			0.25 M
Lysis Buffer	10 mM Tris-HCl pH 7.9	0.75 ml 1 M pH 7.9	per 75 ml → per 25 ml 32 ml
	50 mM KCl	3.75 ml 1 M	→ 1.25 ml 1.5
	1 mM EDTA	0.15 ml 0.5 M	→ 0.05 ml 12 μl
	1 mM PMSF	0.75 ml 100 mM	→ 0.25 ml 0.3
	0.5% Tween 20	1.875 ml 20%	→ 0.625 ml 0.75
	0.5% Nonidet P40	1.875 ml 20%	→ 0.625 ml 0.7
Storage Buffer	50 mM Tris-HCl pH 7.9	50 ml 1 M	per 1 L
	50 mM KCl	50 ml 1 M	
	0.1 mM EDTA	0.2 ml 0.5 M	
	1 mM DTT	1 ml 1 M	
	0.5 mM PMSF	5 ml 100 mM (omit from final)	
	50% glycerol	500 ml	

Grow overnight culture of *E. coli* INV1αF'/p*Taq* in L-broth + ampicillin
 Inoculate 0.25 ml overnight culture into 500 ml L-broth supplemented with 80 mg/L
 ampicillin in 2 L baffle flask at 37°C

Grow ~ 11 hr to an OD₆₀₀ ~0.8 = 0.5 ml (= 0.25 ml of 1 ml IPTG)

Add IPTG to 125 mg/L. Continue growing for 12 hr.
 Harvest cells by centrifugation (8K, 15 min). (2 purple capped 500 ml bottles ~ 250 ml in each.)

Wash in 25 ml **Buffer A**, centrifuge as above.
 Resuspend pellets in 25 ml **Buffer A** + 4mg/ml lysozyme ~ 12.5 ml/pellet

Incubate 15 min at RT
 Add 25 ml **Lysis Buffer**. Transfer to a Pyrex flask. Incubate 75°C 1 hr.

Transfer to centrifuge bottles. Fuge 15,000 rpm, 10 min, 4°C.
 Transfer clarified lysate to a clean Pyrex flask, get stirring with stir bar at RT.

Add 15 gm (NH₄)₂SO₄ (powdered) to stirring lysate at RT.
 Transfer to centrifuge bottles. Fuge 15,000 rpm 10 min.

Harvest precipitate from pellet and surface.
 Resuspend in 10 ml **Buffer A**.

Dialyze against storage buffer for 12 hr at 4°C.
 Change buffer and repeat dialysis step Don't add PMSF to dialysis buffer.
 Dilute 1:1 with sterilized Storage Buffer (No PMSF) and store in small (25 μl)
 aliquots at -80°C.