

PROTOCOL FOR SEQUENCING BY THE MAXAM-GILBERT TECHNIQUE

Below, we describe in abbreviated form the chemical reactions developed by A. Maxam and W. Gilbert for the base-specific modification and cleavage of DNA. A more-detailed description, together with a thorough discussion of the methods used to isolate asymmetrically labeled fragments of DNA and to prepare and run sequencing gels has been published by Maxam and Gilbert in *Methods in Enzymology* (1980) **65** (part 1): 497-559.

Reagents

Dimethylsulfate (DMS) (Aldrich Chemical Co.).

Hydrazine (HZ) (Eastman Kodak).

Formic acid.

Piperidine (Fisher Scientific). Stock solution is 10 M. Dilute to 1.0 M just before use.

95% *Ethanol*.

70% *Ethanol*.

Distilled H₂O.

1 M *Acetic acid*.

0.3 M *Sodium acetate* (pH 5.2).

5 M *NaCl*.

1.2 N *NaOH*.

1 mM *EDTA*.

tRNA. Stock solution is 1 mg/ml in distilled H₂O.

Buffers

DMS buffer

50 mM sodium cacodylate (pH 8.0)
1 mM EDTA

It is usually unnecessary to adjust the pH.

DMS stop

1.5 M sodium acetate (pH 7.0)
1.0 M mercaptoethanol
100 $\mu\text{g/ml}$ tRNA

HZ stop

0.3 M sodium acetate
0.1 mM EDTA
25 $\mu\text{g/ml}$ tRNA

Loading buffer

80% (v/v) deionized or recrystallized formamide
50 mM Tris-borate (pH 8.3)
1 mM EDTA
0.1 (w/v) xylene cyanol
0.1% (w/v) bromophenol blue

Note

Filtration of the above buffers is usually not necessary. However, if a solution is turbid, filter it through a 0.45- μm nitrocellulose filter.

TABLE A.6. SUMMARY OF BASE-SPECIFIC REACTIONS FOR SEQUENCING END-LABELED DNA

	G	G & A	T & C	C	A > C
<i>Mix</i>	200 μ l DMS buffer	10 μ l H ₂ O	10 μ l H ₂ O	15 μ l 5 M NaCl	100 μ l { 2 N NaOH 1 mM EDTA
<i>Chill to</i>	5 μ l [³² P]DNA 0°C	10 μ l [³² P]DNA 0°C	10 μ l [³² P]DNA 0°C	5 μ l [³² P]DNA 0°C	5 μ l [³² P]DNA Heat to 90°C, 6 min.
<i>Add</i>	1 μ l DMS	50 μ l formic acid	30 μ l HZ	30 μ l HZ	150 μ l 1 N acetic acid 5 μ l tRNA (1 mg/ml) 750 μ l 95% ethanol
<i>Incubate</i>	20°C, 3-4 min.	20°C, 5 min.	20°C, 5 min.	20°C, 5 min.	
<i>Add</i>	50 μ l DMS stop 750 μ l ethanol	180 μ l HZ stop 750 μ l ethanol	200 μ l HZ stop 750 μ l ethanol	200 μ l HZ stop 750 μ l ethanol	
<i>Store</i>	-70°C, 10-15 min.	-70°C, 10-15 min.	-70°C, 10-15 min.	-70°C, 10-15 min.	-70°C, 10-15 min.
<i>Centrifuge</i>	10 min.	10 min.	10 min.	10 min.	10 min.
<i>To pellet add</i>	250 μ l 0.3 M NaAc	250 μ l 0.3 M NaAc	250 μ l 0.3 M NaAc	250 μ l 0.3 M NaAc	250 μ l 0.3 M NaAc
<i>Store</i>	750 μ l ethanol -70°C, 10-15 min.	750 μ l ethanol -70°C, 10-15 min.	750 μ l ethanol -70°C, 10-15 min.	750 μ l ethanol -70°C, 10-15 min.	750 μ l ethanol -70°C, 10-15 min.
<i>Centrifuge</i>	10 min.	10 min.	10 min.	10 min.	10 min.
<i>Rinse pellet with</i>	70% ethanol	70% ethanol	70% ethanol	70% ethanol	70% ethanol
<i>Vacuum dry</i>					
<i>To pellet add</i>	100 μ l 1.0 M piper- idine	100 μ l 1.0 M piper- idine	100 μ l 1.0 M piper- idine	100 μ l 1.0 M piper- idine	100 μ l 1.0 M piper- idine
<i>Heat to</i>	90°C, 30 min.	90°C, 30 min.	90°C, 30 min.	90°C, 30 min.	90°C, 30 min.
<i>Lyophilize</i>					
<i>Add</i>	10 μ l H ₂ O	10 μ l H ₂ O	10 μ l H ₂ O	10 μ l H ₂ O	10 μ l H ₂ O
<i>Lyophilize</i>					
<i>Add</i>	10 μ l H ₂ O	10 μ l H ₂ O	10 μ l H ₂ O	10 μ l H ₂ O	10 μ l H ₂ O
<i>Lyophilize</i>					
<i>Add</i>	10 μ l loading buffer	10 μ l loading buffer	10 μ l loading buffer	10 μ l loading buffer	10 μ l loading buffer
<i>Vortex</i>					
<i>Heat to</i>	90°C, 1 min.	90°C, 1 min.	90°C, 1 min.	90°C, 1 min.	90°C, 1 min.
<i>Chill in ice</i>					
<i>Load onto gel</i>					

Reactions should be carried out in siliconized Eppendorf tubes.

Sequencing Gels*20% Acrylamide*

acrylamide	96.5 g
methylene-bis-acrylamide	3.35 g
ultra-pure urea	233.5 g
5× TBE	100 ml
H ₂ O	to 500 ml

Urea mix

urea	233.5 g
5× TBE	100 ml
H ₂ O	to 500 ml

5× TBE

Tris base	54 g
boric acid	27.5 g
0.5 M EDTA (pH 8.0)	20 ml

To make an 8% sequencing gel, mix in a small flask:

20% acrylamide	20 ml
urea mix	30 ml
10% ammonium persulfate (freshly dissolved in water)	0.4 ml

Pour the solution into the barrel of a 50-ml syringe and add 50 μ l of TEMED. Mix rapidly and inject the contents of the syringe (no needle should be used) into a preformed, sequencing-gel mold.