

## RT PCR PROTOCOL

### I. RT/ cDNA synthesis:

Use 1 to 1.5  $\mu\text{g}$  of RNA and add the following reagents to give a total volume of 20  $\mu\text{l}$ :

25mM MgCl <sub>2</sub>	4 $\mu\text{l}$	
10x PCR buffer (PE)	2 $\mu\text{l}$	
DEPC tx H <sub>2</sub> O	1 $\mu\text{l}$	DTT 2 $\mu\text{l}$
10mM dNTP's (mix)	8 $\mu\text{l}$	
RNAase IN	25U	65 $\mu\text{l}$
(*) 3' Primer	200ng	
DEPC tx H <sub>2</sub> O to total volume of 20 $\mu\text{l}$ .		

(\*) 3' primer could be a random primer (hexamer), oligo dT primer, or a specific 3' primer.  
Reverse transcriptase used is provided by PE in their RT PCR Kit.

Assemble above reagents in 0.5ml tube, incubate at :  
42° C 60 mins.; 99° C 5mins.; 4° C 5 mins.

### Amplification/PCR Protocol:

25mM MgCl <sub>2</sub>	4 $\mu\text{l}$
10x PCR buffer (PE)	8 $\mu\text{l}$
5' Primer	200ng
(*) 3' Primer	200ng
Taq Polymerase	2.5 U
DEPC tx H <sub>2</sub> O to total volume of 80 $\mu\text{l}$ .	

(\*) If the 3' primer was used in the 1st reaction it is not necessary to add any more primer.

Cycle = 1	2 mins. 97° C
Cycle = 40	1 min. 95° C; 1 min. 37 to 65° C (check primer)
Cycle = 1	7 mins. 37 to 65° C

Check product on a 1 - 1.5% agarose gel.