

Sequencing

I. Denature DNA - for double-stranded rxns.

To each tube add:

2 μ l DNA ($\sim 1.5 \mu$ g/ μ l)
2 μ l 2M NaOH + trace EDTA
6 μ l dH₂O

- 5 min. @ 37° (not more)

Then add:

60 μ l EtOH (~~100%~~^{95%})
10 μ l 1M NaOAc

- mix EtOH + NaOAc, then add together
- vortex immediately + put in dry ice/EtOH bath for 10'
- spin cold 15'; take off sv
- 70% EtOH rinse
- speed-vac dry

II. Annealing primer.

to each tube containing dry pellet, add:

2 μ l "Rx buffer"
2 μ l primer
6 μ l H₂O

- add H₂O + buffer first
- immediately add primer
- 30' @ 37°
- 10' (up to 4-5 h.) RT
- if using SS DNA, put instead in 65° bath + cool to RT ($\sim 30-45$ min.); same ingredients

III. Reactions

Prepare plate:

2.5 μ l termination mixes to each well (ddA, C, G, T) - for dGTP or dTTP.

Cover + set aside - otherwise will evaporate!

Then add:

1.0 μ l DTT
1.4 μ l H₂O
0.4 μ l label mix
0.7 μ l dATP (³⁵S)
0.25 μ l enzyme
1.75 μ l enz. dilution buffer

- mix DTT, H₂O, label mix (dGTP or dTTP) + ³⁵S dATP; put on ice
- add enzyme separately, to cold dilution buffer, then to rest of stock sol'n. Use ASAP!
- add 5.5 μ l stock to each tube from above; leave in tip
- add 3.5 μ l from each tube into each well on plate

Sequencing, cont.

- ... • cover + spin down to mix contents
- 37° for 3-5'
- add stop sol'n, 5.5 μ l to each well
- cover + freeze
- before loading, heat at least 2-3'.
- total time from adding stock sol'n (5.5 μ l) + putting in 37° bath should be 2-5'.

IV. Gel

5% acrylamide - 1L. (for long runs)	7% (short + long)	6% (... a compromise...)
47.5 g acrylamide	66.5 g	57 g
2.5 g bis	3.5 g	3 g
0.5 x TBE	0.5 x	0.5 x
498 g (480g?) urea (= 8.3M)	498 g	498 g

immed. before pouring:

30 μ l TEMED
500 μ l 10% APS

Polymerize \geq 1 h.

To pre-run: use 0.5x TBE (or other buffer); make sure air is out of lower + upper areas between plates. Run at 35-50 mAmps; \approx 1000 to 2000 + V. Should be warm to the touch; pre-run 45 min - 1 h., or until warm + over \sim 1500 V.

Run; short runs \sim 1h 45 m; long runs \sim 3h 30 min.
Before running, force urea out of top loading area. After loading, rm. bubbles
take down + fix in 10% EtOH / 10% Glacial acetic acid.
Dry, + put on film.