Sample Labeling

Step 1: Fluorescent Labeling of gDNA Kit: Agilent SureTag DNA Labeling Kit

- 1. Set the PCR thermocycler on the AGLB program and allow the block to heat up to 94 degrees Celsius.
- 2. Thaw the samples on ice.
- 3. Dilute all of the test samples and reference DNA in a separate 1.5 ml eppendorf tube to the lowest concentration (20ul) of the samples with nuclease-free water.
- 4. From the diluted samples calculate the volume of DNA to obtain 1 ug of DNA. Add the appropriate volume in an 8 strip PCR tube containing 2.5 ul of random primers and the appropriate amount of water to a total volume of 15.5 ul of total reaction volume per sample. Place a cap and spin the reaction mix down.
- 5. For the reference sample labelling, if you have 8 test samples, then do separate labelling reaction for 8 reference samples. Place a cap on it and spin the samples down.
- 6. Place the samples in the PCR thermocycler and run the incubation at 94 degrees Celsius for 10 minutes and then 4 degrees indefinitely, until ready to proceed to the next step.
- 7. DO STEP 7 IN DARK. While it's incubating, thaw the reagents listed below. Prepare two labeling master mixes, one for the test samples (Cyanine-3 dUTP's) and the other for the reference sample (Cyanine-5 dUTP) in the following manner per sample in separate tubes. For 8 samples, make2 separate master mix for 9 reactions, one using Cy3 and the other using Cy5. Store them on ice until ready.

Reagents	x 1
5X buffer	5 ul
10 X dNTP	2.5 ul
Cy3/Cy5	1.5 ul
Exo(-) Klenow	0.5 ul

- 8. Remove the samples from the thermocycler and spin the samples down. Remove the caps for the test samples and add 9.5 ul of the master mix containing Cyanine-3 dUTP to each of the samples. Place a fresh cap and spin them down.
- 9. Repeat step 8 with the reference samples, except add 9.5 ul of the master mix containing Cyanine-5 dUTP.

Place the tubes in the PCR thermocycler and proceed to the next step in the program. Incubate the samples at 37 degrees Celsius for 2 hours, 65 degrees Celsius for 10 minutes, and 4 degrees Celsius indefinitely.