

Stripping and Reprobing Membranes in Western Blotting

1. Wash Membrane Once in TBST for 5 mins.
2. Submerge the membrane in Stripping Buffer
 - 100mM 2- Mercaptoethanol - 0.7 mls (100% BME)
 - 2% SDS - 20 mls (10% SDS)
 - 62.5mM Tris-HCl pH 6.7 - 3.125 mls (2M Tris pH 6.7)

to 100 mls Strip solution.
3. Incubate at 50°C for 30 mins. For more stringent conditions incubate at 70°C
4. Wash Membrane 2 x for 10 minutes