

Sucrose Gradients

Buffer A: Extraction Buffer

100mM Tris pH 9.0 - 1ml 1M Tris pH 9.0

100mM NaCl - 1ml 1M NaCl

5mM KCl - 50ul 1M KCl

0.5mM MgCl₂ - 5ul 1M MgCl₂

1mM CaCl₂ - 10ul 1M CaCl₂

0.5% NP40 - 500ul 10% NP40

ddH₂O - to 10 ml * Add 100ul PMSF (100mM) just before use

Buffer B: Gradient Buffer

10mM HEPES pH 7.8 - 500ul 1M HEPES pH 7.8

5mM Na₂HPO₄ - 250ul 1M Na₂HPO₄

5mM KCl - 250ul 1M KCl

0.5mM MgCl₂ - 25ul 1M MgCl₂

1mM DTT - 500ul 100mM DTT

ddH₂O - to 10 ml

- Pellet cells 50-100 million cells 1500 rpm for 5 mins
- Wash cells once in 25mls of cold PBS and centrifuge as above
- Aspirate cells to dryness keeping on ice
- Add 200ul of Buffer A and resuspend cells thoroughly. Transfer to douncer and dounce 20-30 strokes on ice.
- Transfer to 1.5ml eppendorf and allow to extract on ice for 20 minutes.
- Centrifuge at 4°C for 20 mins 45,000 xg (13,200 rpm) to remove cell debris and nuclei.
- Layer supernatant gently and immediately onto pre-cooled 5-30% sucrose gradients (12 tubes each tube)
- * - Set up gradient using peristaltic pump. Continuous gradient 5-30%

5% sucrose - 5g/100mls Buffer B

30% sucrose - 30g/100mls Buffer B

- Centrifuge in SW41 swinging bucket rotor at 4°C.

260,000 xg for 6 hours

110,000 xg for 14 hours

- Collect 14 fractions and the pellet

