

T CELL ROSETTING WITH AET/SRBC

AET = 2-amino-ethyl isothiuronium bromide hydrobromide

SRBC = sheep red blood cells

PREPARATION OF AET/SRBC

1. Wash qs for 10mls packed SRBCs with cold, sterile normal saline in a 50ml tube 4 x 10' ea. @ 1500rpm.
2. For each 2mls of packed SRBCs weigh 0.40 grams AET and add to 10mls water. Adjust to pH 8.9 w/ 10N NaOH and sterile filter.
3. Add 2.0mls packed SRBCs to 8.0mls AET solution in a 50.0ml tube.
4. Mix gently and incubate @ 37 C x 15'.
5. Wash 4 x 10' ea. w/ cold saline and once more w/ cold RPMI or balanced salt solution + 10% FCS.
6. Resuspend in approx. 50mls cold R10 and count a 1:1000 dilution.
7. Store AET/SRBCs @ 4 C up to one week.

T CELL ROSETTING

0. ALL STEPS @ 4 DEGREES C.
1. Mix 1 x 10⁸ buffy coat MNCs + 1.5 x 10¹⁰ AET/SRBCs in a 50ml tube. Spin x 5' @ 1000rpm @ 4 C.
2. Incubate @ 4 C x 1.5hrs to o.n.
3. Dilute to 45mls w/ R10. Pipette w/ 25ml pipette to separate rosettes from non-Ts.
4. Layer 15mls over 15mls ficol-hypaque. Spin 20' @ 1200rpm, increasing to 3000rpm for the final 5'.
5. Recover the non-Ts from the interface and wash 3 x 10' ea. w/ R10.*
6. Expect 10-20% recovery of number rosetted.

* Lyse contaminating SRBC's by suspending 15⁻ pellet in 3-5mls Tris-NH₄Cl for 3-5min. Suspension will change from opaque to translucent red. continue w/ ~~up~~ step 5 washes. Pellet should be white. Repeat if needed.

a) Tris : 20.6 gm Tris base / l pH 7.2 - 7.4

b) NH₄Cl: 8.3 gm / l

c) Working Tris/NH₄Cl A:B = 1:9