

<b>Department of Medicine</b> <b>Division of Hematology-Oncology</b> <b>Stem Cell and Xenograft Core</b> <b>Standard Operating Procedure</b>		
<b>SOP No.: ---</b>	<b>Title:</b> <b>Thawing Viably Frozen Cells</b>	Page 1 of 1
Revision No. -		
Effective Date: --/--/----		
Originator: (sign and date)		
Unit Supervisor: (sign and date)		

## 1.0 PURPOSE

Instructions for thawing viably frozen cells

## 2.0 MATERIALS

50 mL conical tubes (polypropylene) - *Catalog# 14-959-49A, Supplier – Fisher*  
15 mL conical tubes (polypropylene) - *Catalog# 14-959-70C, Supplier - Fisher*  
Cell strainer - *catalog# 22-363-546, supplier – Fisher*  
PBS – *Catalog# 14190-235, Supplier – Fisher*  
Fetal Bovine Serum – FBS

## 3.0 PROCEDURE

- 3.1 Prepare 40mL PBS containing 2% FBS (PBS2) in a 50mL conical tube.
- 3.2 Place PBS2 in fridge overnight to chill to 4C.
- 3.3 Take cryotube out of the freezer, spray with 70% ethanol.
- 3.4 Add 1mL cold PBS to top of frozen cells in cryotube.
- 3.5 Wait a few minutes to allow cells on top to thaw.
- 3.6 Keep thawing cells by pipetting back and forth PBS2 between conical tube and cryotube.
- 3.7 Continue to transfer thawed cells to conical tube.
- 3.8 Once all cells transferred to conical, centrifuge at 300xg (1200 RPM) for 5 minutes.
- 3.9 Discard supernatant.
- 3.10 Resuspend pellet in 50mL conical tube to appropriate volume for cell counting.
- 3.11 Filter cell using cell strainer to remove dead clumped cells if necessary.
- 3.12 Counted Cells...see Cell Counting Protocol.
- 3.13 Centrifuge 300xg=1200RPM for 5 minutes.
- 3.14 Discard supernatant.
- 3.15 Resuspend cells in PBS2 and use as needed.