

Robert A. Thomas, 07:28 PM 8/15/00 , Re: Hello

Date: Tue, 15 Aug 2000 19:28:02 -0400 (EDT)
From: "Robert A. Thomas" <rathomas@cmb.biosci.wayne.edu>
To: Erle Robertson <esrobert@umich.edu>
Subject: Re: Hello

Here is the recipe and protocol that you requested:

THE CELL LYSIS BUFFER IS COMPRISED OF TWO BUFFERS:

BASIC CELL LYSIS BUFFER The Basic buffer is made as a 5X stock which is autoclaved

15 mM CaCl ₂ (MW 111.0)	0.167g/100 ml
10 mM Magnesium acetate (MW 214.46)	0.215g/100 ml
0.5 mM EDTA (w/ Na ₂ MW 372.2)	0.019g/100 ml
50 mM TRIS-HCl [pH 8.0] (MW 121.14)	0.606g/100 ml

COMPLETE CELL LYSIS BUFFER is made as follows:

1X BASIC BUFFER (See above)	20 ml /100 ml
1 mM DTT (MW 154.2)	0.154g /100ml
0.5% v/v NP40	0.5 ml/100 ml

CELL LYSIS PROTOCOL

- 1) Transfer cells suspension into 15 centrifuge tube and gently pellet cells.
- 2) Wash cells 2X with 2X DPBS (save supernatants)
- 3) Resuspend the pellets in 125 ul of COMPLETE CELL LYSIS BUFFER.
- 4) Incubate 20' at room temperature with vortexing
- 5) Store cell lysates at -80 until ready for assay (PCR, luciferase, CAT)

If you have any questions, let me know.

Bob

Robert Anthony Thomas	Phone: 313-577-3555 (Off)
Department of Pediatrics	Phone: 313-577-0405 (Lab)
Wayne State University	Fax: 313-577-6200
5047 Gullen Mall, Detroit, MI 48202	E-mail: rathomas@cmb.biosci.wayne.edu

On Tue, 15 Aug 2000, Erle Robertson wrote:

>
> Hi Robert,
> What's going on buddy. I tried calling you a couple of times with no
> response. My student will do the transfection tomorrow in both the 8E5 and
> ACH2 cell lines and we need to know what the protocol for preparing the
> cells for PCR analysis of the HIV transcripts from the LTR. Also we need
> some p24 and Tat antibodies for western blot to correlate these studies
> with the RNA.
>
> Thanks and I look forward to hearing from you again,
> Erle
>
> Erle S. Robertson Ph.D.
> Assistant Professor
> Microbiology and Immunology and
> Comprehensive Cancer Center
> University of Michigan Medical School
> 3217 CCGC (office)
> 3111 CCGC (laboratory)
> Ann Arbor, MI. 48109
> fax: 734-647-9647
> phone: 734-647-7296
> Email: esrobert@umich.edu