

Jae Jung, 05:20 PM 1/7/98 -, Re: HVS

Date: Wed, 7 Jan 1998 17:20:44 -0500  
From: jjung@warren.med.harvard.edu (Jae Jung)  
Subject: Re: HVS  
To: Erle Robertson <esrobert@umich.edu>

*HVS Transformation  
of Human PBMC's*

Dear erle,  
happy new year. have you sent the same mail before? that means  
i did not reply back. sorry about that. life is still hectic here.  
papers, renewal and other stuffs! it is good to hear that things are  
moving now. since we are making il-2 independent T cell transformation,  
protocol we are using is very simple.  
Primary peripheral blood mononuclear cells (PBMCs) from common marmosets  
Callithrix jacchus) were purified from 3 ml heparinized blood specimens  
using lymphocyte separation medium (Organon Teknika Corp., Malvern, Pa.).  
Washed cells were resuspended and cultured in RPMI 1640 medium supplemented  
with 20 % (vol/vol) heat inactivated fetal bovine serum and 5 µg/ml of beta  
-mercaptoethanol. HVS infection of marmoset PBMCs was performed in 12 well  
tissue culture plates by addition of 1ml of HVS. Immortalization  
of primary lymphocytes to IL-2 independent growth was typically established  
within a month after infection. Uninfected common marmoset PBMCs were also  
cultivated by stimulation with 1 µg of phytohemagglutinin per ml for 48 hrs  
followed by culture in RPMI 1640 with 20% fetal bovine serum and 10% IL-2.  
good luck with your experiments.

*⇒ conc. BME 14.3µM  
5µg/ml = 1.1µg/µl  
Add 500µl / 500µl Med  
R20*

best regards  
jj

Reply Separator

Subject: HVS  
Author: Erle Robertson <esrobert@umich.edu> at HMS-Internet  
Date: 1/7/98 11:51 AM

Dear Jae,

Happy New Year. I hope that the research is going well and that you are  
having a wonderful time with your family and friends. I am doing fine and  
the lab is just beginning to look like experiments are getting done the way  
I would like to see things moving.

I am trying to do some transformation assays with HVS which you supplied  
and I am having trouble. I was wondering if you can supply me with a  
protocol you have used in the past for doing these assays. I would be very  
grateful if you would assist in this manner. Fleckenstein group does not  
want to give up the cosmids for HVS just yet. I am in the process of doing  
the cloning myself and I hope that everything goes well. I am presently  
making a vector so that I can take out unique fragments as they were cloned.

I look forward to hearing from you,

Best and sincere regards,  
Erle

Erle S. Robertson Ph.D.  
Assistant Professor  
Dept. of Microbiology and Immunology,  
University of Michigan Medical School  
Medical Science Building II, Room 6724  
Ann Arbor, MI 48109-0620  
Fax: 313(734)-764-3562  
Phone: 313(734)-647-7296  
Email: esrobert@umich.edu  
PS. New area code for 1998