

Immunofluorescence



	Initials	Date
Prepared By	ER	
Approved By		

WILSON JONES COMPANY G7506 ColumnWrite MADE IN U.S.A.

- A Make Slides
1. Aliquot Cells
 - Microfuge scoops for 5 mins at R. temp.
 - Aspirate soup using a pasteur pipette adapted with a pipette tip
 - add 1 ml of 1X PBS and pipette cells up and down to mix (Do NOT vortex)
 - Aspirate leaving approximately 20ul of 1X PBS
 - mix cells by pipetting to uniformly mix cells.
 2. Add approx. 1ul of cells to cells and spread evenly to distribute cells.
 - let cells air dry.
 3. Fix in 1:1 Acetone/Methanol at -20°C for 10 mins
 4. Air Dry
 5. Block Incubate with 20% Goat Serum 30 mins at room temperature. (Humidity Chamber).
 6. Wash with 1X PBS
 7. Incubate (Humidity Chamber) with Specific antibodies
 - 1. α Bam Z
 - 2. α gp110

} 1/200 - 1/500 dilute in PBS 1X

 - incubate 1-2 hrs at room temperature.
 8. Wash well with 1X PBS 4-5 5 minute washes.
 9. Incubate with Goat α Mouse Ig - FITC 1/2000 - 1/2000
 10. Wash well with 1X PBS.
 11. Keep Moist.
 12. Add Anti Fade 100ul and cover slip. Flip slide and press firmly to remove any bubbles trapped



1. Goat Serum
2. α Bam 2 or α IgP110
3. Goat α Mouse Ig - FITC

1. Goat Serum
2. —
3. Goat α Mouse Ig - FITC

e.g. BTAB, BTAB/Bam 2, BTAB/B95-8

BTAB Bam 2

$$\frac{42 + 44}{2} = 43$$

$$43 \times 2 \times 10^4$$

$$= 8.6 \times 10^5$$

$$= 860,000 \text{ cells/ml}$$

40 μ l cells + 40 μ l trypan blue = $\frac{1}{2}$ df

	# 1	# 2	# 3	# 4
Volume	2.0ml	1.5ml	1ml	0.5ml
# cells	1.7mill.	1.3mill	900,000	500,000 / 20 μ l 1x PBS
# cells/ μ l	85,000	65,000	45,000	25,000 / 1 μ l 1x PBS