

VIRUS PREPARATION

1. Remove all 10-20cc from T25 after 3-4 days incubation and spin ~1500rpm (~3.5 on control knob) x 5 min. in 50cc tubes at 4 degrees C.
2. Draw off supernatant with a 10 or 30cc syringe (18g needle) leaving 1-2cc supernatant with cell pellet and place syringe (now with the supernatant) on ice.
3. Take the tube with the cell pellet and 1-2cc of residual RPMI (see #2) and freeze/thaw cell pellet x 3 in dry ice/ethanol bath and 37 degree C. bath (tubes can be conveniently held in test tube rack for freeze/thaw). Vortex tubes after each thaw to mix.
4. Spin 2500 rpm (~5 on knob) x 5 min. at 4 degrees C. and use the syringe already holding the virus supernatant (see #2) to collect the residual 1-2cc supernatant (avoid taking cell fragments)
5. Filter virus with 0.45um filter to get rid of cell fragments.
- (6. (Optional- steps 6-8, especially useful if transfection was incubated in 20cc RPMI rather than 10cc) Spin 8500rpm x 2 hrs at 4 degrees C. in sterile Beckman tubes and mark tube where virus should pellet.)
- (7. Aspirate off supernatant but leave ~1.5-4cc. Let sit on ice in cold room overnight to resuspend virus and in the morning gently pipette up and down to resuspend.)
- (8. In morning, infection can be done using concentrated virus (#7)).
9. Part of virus can be frozen and part can be used for infections.