





Attention: New DNA Sequencing Facility Guidelines for Sample Submission

Effective: September 13, 2004

- Starting on Monday, September 13, 2004, researchers will be asked to submit sequencing samples in **strip conjoined 0.2 ml strip tubes**.



- Please note that **No Samples** will be accepted in **1.5 ml tubes** after September 27, 2004.
- If you happen to forget, you will be asked to transfer your samples to 0.2 mL strip tubes in the facility. (Pipettor, tips and tubes will be provided by the facility)
- The tubes are to be numbered consecutively so that they match the corresponding request sheet. The User name and PI name must also be written on the tubes at least once:



Sample sheet:

	Sample Name. Primer Name	Template Type	Primer Added Yes or No?	Standard primer to be added by Facility	Special Instructions (GC rich, etc)
1	Name of Sample 1				
2	Name of Sample 2				
3	Name of Sample 3				