

The following restriction enzymes do not cut ptTS-Neo:

AccI, Afl II, AgeI, AscI, Asp718 I, Bcl I, BlpI, BsiWI, BsmI, BspEI,

BssHII, Bst1107 I, BstEII, DraIII, Eco47 III, EcoRV, FseI, HpaI, KpnI,
MluI,

NotI, NruI, PacI, PmeI, PshAI, PstI, Sal I, SgfI, SgrAI, SmaI,

SrfI, Sse8387 I, SwaI, XcmI, XhoI, XmaI.