
















# Alternate N95 Respirators During COVID-19



Scan to view on Intranet

If You Have Been Fit Tested On:	You Can Also Use:	Or This	Or This	Or This
 <b>3M 1860 Regular</b>	 <b>3M 8210/8210+</b>	 <b>3M 1870+/1870/9210+</b>	 <b>Kimberly Clark Regular</b>	
<b>Fit Testing for Alternate Required?</b>	<b>NO</b>	<b>YES</b>	<b>YES</b>	
 <b>3M 1870+</b>	 <b>3M 1870</b>	 <b>3M 9210+</b>	 <b>3M 1860 Regular</b>	 <b>Kimberly Clark Regular</b>
<b>Fit Testing for Alternate Required?</b>	<b>NO</b>	<b>NO</b>	<b>YES</b>	<b>YES</b>
 <b>3M 1860S</b>	 <b>3M 8110S</b>	 <b>3M 1870/1870+/9210+</b>	 <b>Kimberly Clark Small</b>	 <b>3M 1860 Regular</b>  <b>3M 8210/8210+</b>
<b>Fit Testing for Alternate Required?</b>	<b>NO</b>	<b>YES</b>	<b>YES</b>	<b>YES</b>