

LSR A**SN: H1095****Blue Laser (488nm)**

515/20 Blue **695/40 Blue**
Blue PMT B **Blue PMT A**

AF488	7-AAD
BB515	BB700
CFSE	PerCP-Cy5.5
FITC	PerCP-eF710
GFP	
YFP	

Red Laser (640nm)

660/20 Red **713/54 Red** **780/60 Red**
Red PMT C **Red PMT B** **Red PMT A**

APC	AF680	APC-Cy7
AF647	APC-R700	APC-H7
eF660	Draq5	APC-AF750
eF670	Draq7	APC-eF780
DyLight 649	AF700	APC-FIRE
TO-PRO-3		

Green Laser (532nm)

575/26 Green **610/20 Green** **660/20 Green** **695/40 Green** **780/60 Green**
Green PMT E **Green PMT D** **Green PMT C** **Green PMT B** **Green PMT A**

BYG584	PE-CF594	PE-Cy5	PE-Cy5.5	PE-Cy7
DsRed	PE-Texas Red		PE-AF700	PE-Vio770
PE	PE-AF610			
tdTomato	PE-eF610			
RFP	PI			
RY586	mCherry			

Violet Laser (405nm)

450/50 Violet **515/20 Violet** **560/40 Violet** **585/42 Violet** **605/40 Violet** **660/40 Violet** **710/50 Violet** **790/60 Violet**
Violet PMT H **Violet PMT G** **Violet PMT F** **Violet PMT E** **Violet PMT D** **Violet PMT C** **Violet PMT B** **Violet PMT A**

AF405	BV510	BV570	BV570	BV605	BV650	BV711	BV785
BV421	Fix L/D Aqua						BV786
SB436	V500						
CellTrace Violet							
DAPI							
eF450							
Pac Blue							
V450							
Sytox Blue							

LEGEND

AF - Alexa Fluor
BB - Brilliant Blue
BUV - Brilliant Ultraviolet
BV - Brilliant Violet
BYG - Brilliant Yellow Green
eF - eFluor
L/D - Fixable Live/Dead