

LSRFortessa CJ2

SN: H64717700062

Blue Laser (488nm)

515/20 Blue 710/50 Blue
Blue PMT B Blue PMT A

AF488	7-AAD
BB515	BB700
CFSE	PerCP-Cy5.5
FITC	PerCP-eF710
GFP	
YFP	

Red Laser (640nm)

670/30 Red 710/50 Red 780/60 Red
Red PMT C Red PMT B Red PMT A

APC	AF680	APC-Cy7
AF647	APC-R700	APC-H7
eF660	Draq5	APC-AF750
eF670	Draq7	APC-eF780
DyLight 649	AF700	APC-FIRE
TO-PRO-3		

Green Laser (532nm)

575/26 Green 610/20 Green 660/40 Green 710/50 Green 780/40 Green
Green PMT E Green PMT D Green PMT C Green PMT B Green PMT A

BYG584	PE-CF594	PE-Cy5	PE-Cy5.5	PE-Cy7
DsRed	PE-Texas Red		PE-AF700	PE-Vio770
PE	PE-AF610			
tdTomato	PE-eF610			
RFP	PI			
RY586	mCherry			

Violet Laser (405nm)

450/50 Violet 525/50 Violet 570/20 Violet 605/30 Violet 655/40 Violet 710/40 Violet
Violet PMT F Violet PMT E Violet PMT D Violet PMT C Violet PMT B Violet PMT A

AF405	BV510	BV570	BV605	BV650	BV711
BV421	Fix L/D Aqua				
SB436	V500				
CellTrace Violet					
DAPI					
eF450					
Pac Blue					
V450					
Sytox Blue					

LEGEND

AF - Alexa Fluor
BB - Brilliant Blue
BUV - Brilliant Ultraviolet
BV - Brilliant Violet
BYG - Brilliant Yellow Green
eF - eFluor
L/D - Fixable Live/Dead